# INTRASPINE®

## DEVICE FOR THE INTER SPINOUS SPACE WITH LAMINAR SUPPORT

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notice d'instructions</a>	Page	6
de	<a href="#">Gebrauchsanweisung</a>	Seite	11
it	<a href="#">Istruzioni per l'uso</a>	Pagina	16
es	<a href="#">Instrucciones de uso</a>	Página	21
pt	<a href="#">Instruções de uso</a>	Página	26
ελ	<a href="#">Οδηγίες χρήσης</a>	σελίδα	31
pl	<a href="#">Instrukcja obsługi</a>	Strona	36
cs	<a href="#">Návod k použití</a>	Strana	41
sl	<a href="#">Navodila za uporabo</a>	Stran	46
ru	<a href="#">Инструкция по применению</a>	Страница	51
kr	<a href="#">사용 지침</a>	페이지	57
en	<a href="#">Patient information leaflet (Australian requirements)</a>	Page	60

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This release is the last update of the instructions of use and replace the previous edition

## IntraSPINE®

### DEVICE FOR THE INTER-SPINOUS SPACE WITH LAMINAR SUPPORT

#### SINGLE USE STERILE PRODUCT

#### **DESCRIPTION**

The IntraSPINE® prosthesis is an interlaminar medical device that comprises a dimethyl siloxane (silicon) wedge (available in several sizes) adapted to the intervertebral space and a ligament in polyethylene terephthalate.

The wedge is covered with woven polyethylene terephthalate with an area covered by sheeting on the front face in the medullary zone.

Tension is applied to the ligament by a titanium loop (TA6V).

#### **MATERIALS**

IntraSPINE® prosthesis: Dimethyl siloxane, Polyethylene terephthalate, Titanium.

IntraSPINE® ancillaries and trial prostheses: Stainless steel, Acetal copolymer, Dimethyl siloxane.

Non-human and non-animal origin – Non-resorbable.

#### **INDICATIONS**

The IntraSPINE® implant is recommended for use in the lumbar spine including in the lumbosacral area (from L1 to S1) in the following cases:

- Arthropathic facet-syndrome
- Foraminal stenosis
- Degenerative discopathy
- Intervertebral ligament insufficiency

#### **PERFORMANCE**

The IntraSPINE® prosthesis is an elastic support device with a cushioning effect. Its interlaminar implantation enables it to be close to the centre of rotation, thereby producing the cushioning effect and providing tension on the interspinous ligamentary structures due to the elasticity of the supra-spinous wedge.

The principal function of the ligaments surrounding the upper and lower vertebrae in relation to the area receiving the prosthesis is to assist the posterior ligament during flexion.

IntraSPINE® has been validated for a long term implantation in the human body.

#### **CONTRAINDICATIONS**

1) Do not use the prosthesis in the following cases:

- Allergy to one of the components
- Growing child
- Infected site
- Pregnant woman

2) Factors likely to compromise successful implantation:

- Severe osteoporosis
- Major deformations of the spinal column
- Local bone tumors
- Systemic or metabolic disorders
- Infectious diseases
- Obesity
- Drug addictions
- Intense physical activity, for example competitive sports or hard labor
- Surgery at more than 2 levels with IntraSPINE Prosthesis

#### **UNDESIRABLE SECONDARY EFFECTS**

**1) Complications likely related to the device:**

- Infection of the prosthesis
- Device breakage (rupture of the ligament, fixed wing breakage, knitting tearing)
- Uncrimping of the needle
- Device loosening permanent elongation of the ligament
- Device migration displacement of the wedge
- Spondylolisthesis to the dynamically stabilized segment
- Persistent or worsening pain attributed to the dynamically stabilized level and required revision surgery with interbody fusion at that level
- Spinous process fractures intraoperative post-operative or during device removal (removal due to infection) and can be associated with dislocation of the implant
- Late reoperation due to absence of recovery
- Disability
- Discomfort
- Wrong level implantation, inappropriate or inaccurate positioning (posterior or wrong level)
- Seroma/swelling around the device
- Inflammatory reaction, allergy, foreign body reaction
- In certain cases the medical instrument had to be removed and pedicle screw fusion had to be performed.

**2) Complications unlikely to be related to the device, but related to surgery or to the patient's conditions:**

- Pseudomeningocele, fistula, dural tear, persistent CSF leakage, meningitis
- Loss of neurological function,
- Cauda equina syndrome, neuropathy, neurological deficits, arachnoiditis, and/or muscle loss, compression around nerves and/or pain
- Urinary retention or loss of bladder control or other types of urological system compromise
- Fracture, Bone loss or decrease in bone density
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: Adjacent segment degeneration
- Non-union or pseudarthrosis, delayed union. Mal-union,
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function
- Hemorrhage, hematoma, edema, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, wound infection, damage to blood vessels, or other types of cardiovascular system compromise
- Development of respiratory problems

#### **MRI SAFETY:**

The implants are composed of non-ferromagnetic materials and present a geometry not likely to generate induced current. Moreover, as they are fixed to tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms

of heating and migration of the implant, has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded a conditional compatibility for MRI scan up to 1.5 Tesla. As a precautionary measure, it is recommended MRI scans be avoided within 48 hours of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential. It should be noted that devices presenting a high contrast with the biological environment can generate "artifacts" that must be taken into account for the perfect execution and interpretation of imaging exams. For this purpose, the patient who has this implant must be advised to warn the health professionals concerned (radiologists and radiology operators) in as far as possible about the presence of this implant before these exams. The Safety of the MR environment has not been tested.

**POSTOPERATIVE PRECAUTIONS**

The patient must be informed of the postoperative cares and of precautions. The use of a brace could be necessary but is not mandatory. This decision is the surgeon one.

**PRECAUTIONS FOR USE**

Before usage, verify the integrity of the prosthesis, trial prostheses, ancillaries and packaging. Do not use in the event of deterioration of the prosthesis and/or the labels, and/or the trial prostheses and/or the ancillaries and/or the packaging. Do not use if the prosthesis has exceeded its use-by date. Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility. When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the prosthesis are equal to the data on the label of the inner pouch. Otherwise the use of the prosthesis is strictly prohibited. Prevent the prosthesis, the trial prostheses and ancillaries from coming in to contact with objects that may alter their surface. Operation conditions: Temperature from +32°C up to 42°C. Implants should be resistant to biologic fluids and the discharge of body tissues, which they contact during use. Therefore, the IntraSPINE® prosthesis has been validated for a long term implantation and, currently, COUSIN BIOTECH has 10 years of feedback concerning these products. The IntraSPINE® prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of IntraSPINE®. The information contained in this notice is necessary but not sufficient to control the surgical technique. The IntraSPINE® is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials. The prescription of the device is decided on by the surgeon, the only person qualified to do so. Topping off is not contra-indicated but special attention is required. IntraSPINE® is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials.

Use of the system of COUSIN BIOTECH ancillaries and specific trial prostheses is required for implantation of the IntraSPINE® prosthesis.


Implants must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature. Transportation: Transportation of the devices should be done in the closed vehicles in accordance with the rules for those transport vehicles. There are no special transportation requirements for the device.

**IMPORTANT : DO NOT REUSE AND DO NOT RESTERILISE THE INTRASPINE® PROSTHESIS**

As indicated on the product labelling, the prosthesis is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

**WARNING ABOUT ANCILLARIES AND TRIAL PROSTHESES**

The IntraSPINE® ancillaries are Class I medical devices, intended for temporary use and are re-usable. The IntraSPINE® trial prostheses are Class IIa medical devices, intended for temporary use and are re-usable. Unlike the IntraSPINE® prosthesis, ancillaries and trial prostheses (included with the instrument kit/set) are furnished unsterile.

**KITACINTRA / RCBININTRU**

**Before first use**

- Ancillaries and trial prostheses must be:
- Cleaned in the appropriate way
  - Sterilised in an autoclave

**DECONTAMINATION AND CLEANING OF ANCILLARIES**

Decontamination and cleaning are performed under the responsibility of the healthcare centers. The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments

**The detachable instruments have to be dismantled**

**Products advised for the cleaning**

- an enzymatic detergent neutral
- Or products adapted for the cleaning (neutral or soft alkaline) of chirurgical instruments in compliance with the regulations in effect.

- **Warning:** Don't use corrosive factor or caustic cleaning product  
**2 process are possible**

**A- Chart Automated cleaning in a disinfecter washer with manual Pre-cleaning:**

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush <b>(Do not utilize metal cleaning brushes)</b>	Until all visible soils removed
	2	Pre-soak	Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.	- Detergent - Water jet pistol	Until product is visually clean
	3	Ultrasonic cleaning	Clean the instrument in the ultrasonic tray including the solution. <b>Be careful every unclench connector or screws and bolts by the vibrations will be reassure.</b>	- Ultrasonic tray	Minimum duration 10 minutes, this duration can vary according product
	4	Soak	Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position.	- Reserve osmosis water	Minimum 30 seconds of soak
	5	Automated Washer	Place entire disassembled device into the automated washer	- Automated Washer - Disassembled Device	Minimum total cycle time: 39 minutes

Instructions of cleaning in washer / disinfecter			
Cycles parameters	Time	Minimum temperature	Type of detergent/water
Pre-cleaning	2 minutes	<45°C / <113°F	Tap water
Cleaning	5 minutes	55°C / 131°F	Detergent
Draining	2 minutes	<45°C / <113°F	Tap water
Thermic draining	5 minutes	90°C / 194°F	Reserve osmosis water
Drying	25 minutes	> 70°C / > 158°F	Not applicable

Finishing Steps	6	Final Rinse	<b>Rinse carefully with distilled water</b>	- Distilled water	Minimum 1 minute
	7	Final Drying	<b>Dry</b> devise utilizing medical quality filtered air	- Medical quality filtered air	Until product is visually dry
	8	Tidying	The instruments have to be placed <b>in the right position</b> inside the trays as described <b>on the position map fixed</b> on the tray for the sterilization.		

**B- Chart Just manual cleaning :**

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush <b>(Do not utilize metal cleaning brushes)</b>	Until all visible soils removed
	2	Drying	Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available	- Non-schedding wipe - Medical quality filtered compressed air	Until product is visually dry
Disinfection Steps	3	Disinfection Application	Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument.	- Detergent	Minimum duration 1 minute, this duration can vary according product
	4	Drying	After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water.	- Disposable cloth - Detergent	Until product is visually dry

Disinfection Steps	5	Manual Disinfection	Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution.	<ul style="list-style-type: none"> <li>- Soft bristle brush</li> <li>- Detergent</li> <li>- Syringe</li> </ul>	Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned
	6	Rinse and finition	Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water	<ul style="list-style-type: none"> <li>- Reserve osmosis water</li> <li>- Syringe</li> </ul>	Minimum 1 minute
	7	Final Rinse	<b>Rinse carefully with distilled water</b>	<ul style="list-style-type: none"> <li>- Distilled water</li> </ul>	Minimum 1 minute
	8	Final Drying	<b>Dry</b> devise utilizing medical quality filtered air	<ul style="list-style-type: none"> <li>- Medical quality filtered air</li> </ul>	Until product is visually dry
	9	Tidying	The instruments have to be placed <b>in the right position</b> inside the trays as described <b>on the position map fixed</b> on the tray for the sterilization.		

**Note:** In case of **suspicion of Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers  
In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

### **STERILIZATION OF ANCILLAIRES**

Sterilization is performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

1- The sterilization in autoclave is to make in a **specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

**NB: It isn't recommended** to realize the sterilization with peelables pouches **autoclavable** without a container  Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

Steam sterilization parameters	Values	Values
Cycle Type	Pre-vacuum	Pre-vacuum
Set Point Temperature	134°C / 273°F	132°C / 269.6°F
Cycle Time	Minimum 3 minutes	Minimum 4 minutes
Dry Time	Minimum 20 minutes	Minimum 20 minutes
Cool Time (inside and outside the autoclave chambers)	Minimum 40 minutes	Minimum 40 minutes

3 -Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrality and the correct operation, of the locking systems (clipping, prehension) without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent sterilization process for which it will be fully responsible.

### **STORAGE PRECAUTIONS**

Ancillaries and trial prostheses must be stored in their dedicated container or equivalent packaging to prevent any deterioration.

Keep in a dry place away from sunlight, at room temperature, in its original packaging.

### **EXPLANTATION AND DISPOSAL OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### **INFORMATION REQUESTS AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### **IMPLANT CARD**

**Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.**

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## Patient information leaflet IntraSPINE®

Item	Patient information leaflet
1	<p>a. <b>The name of the device :</b> IntraSPINE®</p> <p>b. <b>The model of the device :</b> INTRACAL08, INTRACAL10, INTRACAL12, INTRACAL14, INTRACAL16</p>
2	<p>a. <b>The intended purpose of the device :</b>            The IntraSPINE® prosthesis is recommended in the lumbar spine including in the lumbosacral area (from L1 to S1 levels) in the following cases:</p> <ul style="list-style-type: none"> <li>- Arthropathic facet-syndrome</li> <li>- Foraminal stenosis</li> <li>- Degenerative discopathy</li> <li>- Intervertebral ligament insufficiency</li> </ul> <p>b. <b>The kind of patient on whom the device is intended to be used :</b>            The target population is skeletally mature patients.</p>
3	<p><b>Special operating instructions for the use of the device :</b> The IntraSPINE® prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of IntraSPINE®. The information contained in this notice is necessary but not sufficient to control the surgical technique.</p>
4	<p>a. <b>The intended performance of the device :</b>            The IntraSPINE® prosthesis is an elastic support device with a cushioning effect. Its interlaminar implantation enables it to be close to the centre of rotation, thereby producing the cushioning effect and providing tension on the interspinous ligamentary structures due to the elasticity of the supra-spinous wedge.            The principal function of the ligaments surrounding the upper and lower vertebrae in relation to the area receiving the prosthesis is to assist the posterior ligament during flexion.            IntraSPINE® has been validated for a long term implantation in the human body.</p> <p>b. <b>Undesirable side effects :</b></p> <ul style="list-style-type: none"> <li>• <b>Complications likely related to the device</b> <ul style="list-style-type: none"> <li>- Infection of the prosthesis</li> <li>- Device breakage (rupture of the ligament, fixed wing breakage, knitting tearing)</li> <li>- Uncrimping of the needle</li> <li>- Device loosening permanent elongation of the ligament</li> <li>- Device migration displacement of the wedge</li> <li>- Spondylolisthesis to the dynamically stabilized segment</li> <li>- Persistent or worsening pain attributed to the dynamically stabilized level and required revision surgery with interbody fusion at that level</li> <li>- Spinous process fractures intraoperative post-operative or during device removal (removal due to infection) and can be associated with dislocation of the implant</li> <li>- Late reoperation due to absence of recovery</li> </ul> </li> </ul>

This leaflet was revised in March 2023

Page

63/72

	<ul style="list-style-type: none"> <li>- Disability</li> <li>- Discomfort</li> <li>- Wrong level implantation, inappropriate or inaccurate positioning (posterior or wrong level)</li> <li>- Seroma/swelling around the device</li> <li>- Inflammatory reaction, allergy, foreign body reaction</li> <li>- In certain cases the medical instrument had to be removed and pedicle screw fusion had to be performed.</li> <li>• <b>Complications unlikely to be related to the device, but related to surgery or to the patient's conditions</b></li> <li>- Pseudomeningocele, fistula, dural tear, persistent CSF leakage, meningitis</li> <li>- Loss of neurological function,</li> <li>- Cauda equina syndrome, neuropathy, neurological deficits, arachnoiditis, and/or muscle loss, compression around nerves and/or pain</li> <li>- Urinary retention or loss of bladder control or other types of urological system compromise</li> <li>- Fracture, Bone loss or decrease in bone density</li> <li>- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: Adjacent segment degeneration</li> <li>- Non-union or pseudarthrosis, delayed union. Mal-union,</li> <li>- Cessation of any potential growth of the operated portion of the spine</li> <li>- Loss of or increase in spinal mobility or function</li> <li>- Hemorrhage, hematoma, edema, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, wound infection, damage to blood vessels, or other types of cardiovascular system compromise</li> <li>- Development of respiratory problems</li> </ul>
5	<p><b>When IntraSPINE® Should Not be Used :</b></p> <p>1) Do not use the prosthesis in the following cases: Allergy to one of the components, Growing child, Infected site, Pregnant woman</p> <p>2) Factors likely to compromise successful implantation: Severe osteoporosis, Major deformations of the spinal column, Local bone tumors, Systemic or metabolic disorders, Infectious diseases, Obesity, Drug addictions, Intense physical activity, for example competitive sports or hard labor, Surgery at more than 2 levels with IntraSPINE® Prosthesis</p>
6	<p><b>MRI Information :</b> The implants are composed of non-ferromagnetic materials and present a geometry not likely to generate induced current. Moreover, as they are fixed to tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of the implant, has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded a conditional compatibility for MRI scan up to 1.5 Tesla.</p> <p>As a precautionary measure, it is recommended MRI scans be avoided within 48 hours of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.</p> <p>It should be noted that devices presenting a high contrast with the biological environment can generate "artifacts" that must be taken into account for the perfect execution and interpretation of imaging exams.</p>

6	<p>For this purpose, the patient who has this implant must be advised to warn the health professionals concerned (radiologists and radiology operators) in as far as possible about the presence of this implant before these exams. The Safety of the MR environment has not been tested.</p>
7	<p><b>Precautions :</b></p> <p>Before usage, verify the integrity of the prosthesis, trial prostheses, ancillaries and packaging. Do not use in the event of deterioration of the prosthesis and/or the labels, and/or the trial prostheses and/or the ancillaries and/or the packaging. Do not use if the prosthesis has exceeded its use-by date. Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility. When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the prosthesis are equal to the data on the label of the inner pouch. Otherwise the use of the prosthesis is strictly prohibited. Prevent the prosthesis, the trial prostheses and ancillaries from coming in to contact with objects that may alter their surface. Operation conditions: Temperature from +32°C up to 42°C. Implants should be resistant to biologic fluids and the discharge of body tissues, which they contact during use. Therefore, the IntraSPINE® prosthesis has been validated for a long term implantation and, currently, COUSIN BIOTECH has 10 years of feedback concerning these products. The IntraSPINE® prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of IntraSPINE®. The information contained in this notice is necessary but not sufficient to control the surgical technique. The IntraSPINE® is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials. The prescription of the device is decided on by the surgeon, the only person qualified to do so. Topping off is not contra-indicated but special attention is required. IntraSPINE® is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials. Use of the system of COUSIN BIOTECH ancillaries and specific trial prostheses is required for implantation of the IntraSPINE® prosthesis. Implants must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature. Transportation: Transportation of the devices should be done in the closed vehicles in accordance with the rules for those transport vehicles. There are no special transportation requirements for the device.</p>
8	<p><b>Materials included in the device:</b> Dimethyl siloxane, Polyethylene terephthalate, Titanium.</p>

This leaflet was revised in March 2023

Page

65/72

Headquarters & Factory: **COUSIN BIOTECH** Allée des Roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –

**Fax** +33 (0) 3 20 14 40 13

[www.cousin-surgery.com](http://www.cousin-surgery.com) - **e-mail** [f.pelletier@cousin-surgery.com](mailto:f.pelletier@cousin-surgery.com)


9	<p><b>Reporting a serious incident :</b>          If you encounter any serious incident with your IntraSPINE® implant you should report it immediately to your doctor. You should also report the incident to COUSIN BIOTECH via the link <a href="mailto:serviceclients@cousin-surgery.com">serviceclients@cousin-surgery.com</a> or <a href="mailto:materiovigilance@cousin-biotech.com">materiovigilance@cousin-biotech.com</a> and to the Therapeutic Goods Administration using the following website: <a href="http://www.tga.gov">http://www.tga.gov</a></p>
10	<p><b>Legal Manufacturer &amp; Manufacturing Facility :</b> COUSIN BIOTECH, Allée des Roses, 59117 Wervicq-Sud, France</p>


en fr de it es pt ελ cs pl sl ru kr	- Symbols used on labelling - Symboles utilisés sur l'étiquette - Bei Etiketten verwandete Symbole - Simboli utilizzati sull'etichetta - Símbolos utilizados en el etiquetaje - Símbolos usados na etiqueta - Σύμβολα που χρησιμοποιούνται στις ετικέτες - Symboly použité na štítku - Symbole użyte na etykiecie - Simboli, ki se uporabljajo pri označevanju - Символы, используемые на маркировке - 라벨 표기사항	 <p>www.cousin-biotech.com/ifu</p>  <p>ENOTICE130 INTRASPINE</p> <p>IFU can be downloaded through the QR code. A hard copy can be sent within 7 days on request by email to ifurequest@cousin-biotech.com or by using the order form on our website.</p>
--	---	---

	en fr de it es pt ελ cs pl sl ru kr	Batch number Numéro de lot Chargennummer Numero di lotto Número de lote Αριθμός παρτίδας Číslo šarže Numer partii Serijska številka Номер партии 제조번호		en fr de it es pt ελ cs pl sl ru kr	Caution (See instructions for use) Attention, voir notice d'instructions Vorsicht (Siehe Gebrauchsanweisung) Attenzione (Vedere manuale istruzioni) Precaución (Veanse las instrucciones de uso ) Atenção, consultar nota de instruções Δείτε τις οδηγίες χρήσης Viz pokyny k použití Pozor: Instrukcja obsługi Pozor (Glejte navodila za uporabo) Смотрите инструкцию по применению 주의(IFU 참조)
	en fr de it es pt ελ cs pl sl ru kr	Reference on the brochure Référence du catalogue Verweis auf der Broschüre Referimento del catalogo Referencia al folleto Referência do catálogo Αριθμός καταλόγου Odkaz v brožurě Nr referencyjny na broszurze Sklic na brošuro Ссылка на брошюру 모델명		en fr de it es pt ελ cs pl sl ru kr	Manufacturer Fabricant Hersteller Produttore Fabricante Fabricante Κατασκευαστής Výrobce Producent Proizvajalec Производитель 제조사
	en fr de it es pt ελ cs pl sl ru kr	Do not reuse Ne pas réutiliser Nicht wiederverwenden Non riutilizzare No volver a utilizar Não reutilizar Να μην επαναχρησιμοποιείται Nepoužívat znovu Nie nadaje się do ponownego użytku Ne uporabljajte več Не предназначен для повторного использования 재사용금지		en fr de it es pt ελ cs pl sl ru kr	Use before: year and month Utiliser jusque : année et mois Verbrauchen bis : Jahr und Monat Utilizzare entro e non oltre: anno e mese Usese antes de: año y mes Utilizar até : ano e mês Ημερομηνία λήξης: έτος και μήνας Spotřebať do: rok a měsíc Należy użyć przed upływem: rok i miesiąc Uporabite pred: leto in mesec Срок годности: год и месяц 사용기한: 년월
	en fr de it es pt ελ cs pl sl ru kr	Keep in a dry place Conserver au sec Trocken lagern Conservare in un luogo asciutto Manténgase en un lugar seco Manter em lugar seco Φυλάξτε το προϊόν μακριά από την υγρασία Uchovujte mimo slunečnmu záření Przechowywać w suchym miejscu Hraniti na suhem mestu Хранить в сухом месте 건조한 곳에 보관		en fr de it es pt ελ cs pl sl ru kr	Do not use if the packaging is damaged Ne pas utiliser si l'emballage est endommagé Nicht verwenden, wenn Packung beschädigt ist Non utilizzare se la confezione è danneggiata No usar si el embalaje está dañado Não usar se a embalagem estiver danificada Να μην χρησιμοποιείται, εάν η συσκευασία έχει υποστεί κάποια βλάβη Nepoužívejte je-li obal poškozený Nie używać jeśli opakowanie jest uszkodzone Ne uporabljajte, če je embalaža poškodovana Не использовать в случае повреждения упаковки 포장 손상시 사용 금지
	en fr de it es pt ελ cs pl sl ru kr	Keep away from sunlight A stocker à l'abri de la lumière du soleil Vor Sonnenlicht geschützt lagern Tenere lontano dalla luce solare Mantérese alejado de la luz solar Manter ao abrigo da luz solar Προφυλάξτε το προϊόν από το ηλιακό φως Uchovujte mimo dosah světla Nie wystawiać na bezpośrednie działanie promieni słonecznych Hraniti ločeno od sončne svetlobe Хранить в темном месте вдали от солнечных лучей 직사광선을 피해서 보관		en fr de it es pt ελ cs pl sl ru kr	Do not re-sterilize Ne pas restériliser Nicht erneut sterilisieren Non risterrilizzare No volver e esterilizar Não reesterilizar Να μην επαναποστειρώνεται Znovu nesterilizujte Nie poddawać ponownej sterylizacji Ne ponovno sterilizirajte Не подлежит повторной стерилизации 재멸균 금지
	en fr de it es pt ελ cs pl sl ru kr	Double sterile barrier system Double barrière stérile Doppelsterilbarriere Doppia barriera sterile Doble barrera estéril Dupla barreira estéril Διπλό αποστειρωμένο σύστημα φραγμού Dvojité sterilní bariéra Podwójna bariera sterylna Dvojina sterilna pregrada двойной стерильный барьер 이중 멸균 장벽 시스템		en fr de it es pt ελ cs pl sl ru kr	Medical device Dispositif médical Medizinisches Gerät Dispositivo medico Dispositivo médico Dispositivo médico Ιατρική συσκευή Zdravotnický prostředek Wyrób medyczny Medicinski pripomoček Медицинское оборудование 의료 기기

	en fr de it es pt ελ cs pl sl ru kr	Manufacturing date Date de fabrication Datum der Herstellung Data di fabbricazione Fecha de fabricación Data de fabricaçã Ημερομηνία κατασκευής Datum výroby Data produkcji Datum izdelave Дата изготовления 제조일자		en fr de it es pt ελ cs pl sl ru kr	Patient name Nom du patient Patientenname Nome del paziente Nombre del paciente Nome do paciente Όνομα ασθενή Jméno pacienta Nazwa pacjenta Ime bolnika Имя пациента 환자 이름
	en fr de it es pt ελ cs pl sl ru kr	Name of practitioner Nom du médecin Name des Praktikers Nome del praticante Nombre del profesional Nome do médico Όνομα ασκούμενου Jméno odborníka Nazwa podmiotu prowadzącego praktykę Ime izvajalca Имя практикующего врача 개업의 이름	<b>UDI-DI</b>	en fr de it es pt ελ cs pl sl ru kr	Unique device identification - device identifier Identification unique du dispositif - identifiant du dispositif Eindeutige Produktidentifikation – Produktkennung Identificazione unica del dispositivo - identificativo del dispositivo Identificación única del dispositivo - identificador del dispositivo Identificação única do dispositivo - número de identificação do dispositivo Μοναδική ταυτοποίηση της συσκευής - αναγνωριστικό της συσκευής Jedinečná identifikace prostředku - identifikační číslo prostředku Unikalna identyfikacja urządzenia - identyfikator urządzenia Edinstveni identifikator pripomočka - identifikator pripomočka Уникальный идентификатор изделия - Идентификатор изделия 고유 기기 식별 기기 식별자
	en fr de it es pt ελ cs pl sl ru kr	Date of surgical procedure Date de l'intervention Datum des chirurgischen Eingriffs Data della procedura chirurgica Fecha de la intervención quirúrgica Data do procedimento cirúrgico Ημερομηνία χειρουργικής επέμβασης Datum kirurškého zákroku Data zabiegu chirurgicznego Datum kirurškega posega Дата хирургического вмешательства 수술 날짜		en fr de it es pt ελ cs pl sl ru kr	MR conditional IRM sous conditions Bedingt MRT-kompatibel MRI con riserva IRM en condiciones IRM sob condições MRI υπό προϋποθέσεις MRI za podmínek Warunkowo bezpieczny w środowisku RM Slikanje z magnetno resonanco pod določeni pogoji Условно совместимо с МРТ MRI-조건부
	en fr de it es pt ελ cs pl sl ru kr	Unique device identification Numéro d'identifiant unique Einmalige Identifikationsnummer Identificazione unica del dispositivo Número de identificación única Número de identificação único Μοναδικός αριθμός αναγνώρισης Jedinečné identifikační číslo Unikalny numer identyfikacyjny Edinstveni identifikator Уникальный идентификационный номер 고유 식별 번호		en fr de it es pt ελ cs pl sl ru kr	Number of implants in the packaging = Qty: Nombre d'implant dans l'emballage = Qté: Anzahl der Implantate in der Packung = Anz.: Numero di impianti nella confezione = Qtà: Número de implante en el embalaje = Cant.: Número de implantes na embalagem = Qtd: Αριθμός εμφυτευμάτων στην συσκευασία = Ποσότητα: Počet implantátů v balení = Množ.: Liczba implantów w opakowaniu = Ilość: Število vsadkov v pakiranju = Količina: Количество имплантов в упаковке = Кол-во: 포장된 이식물 수 = 수량
	en fr de it es pt ελ cs pl sl ru kr	Not sterile Non stérile Nicht steril Non sterile No estéril Não esterilizado Μη αποστειρωμένο Nesterilní Niesterylno Nesterilno Нестерильно 비 살균성		en fr de it es pt ελ cs pl sl ru kr	Reusable Réutilisable Wiederverwendbar Riutilizzabile Reutilizable Reutilizável Επαναχρησιμοποιήσιμο Znovu Použitelné Do wielokrotnego użytku Za ponovno uporabo Нестерильно 재사용 가능
	en fr de it es pt ελ cs pl sl ru kr	Sterilized by gamma radiation. Produit stérile. Méthode de stérilisation : irradiation Sterilisierungsmethode : Gammastrahlung. Il prodotto è sterile. Sterilizzazione : raggi Gamma. Esterilizado por rayos gamma. Produto estéril. Método de esterilização: irradiação Στείρο προϊόν. Μέθοδος αποστείρωσης: ακτινοβολία Sterilizováno gamma zářením. Steryliżowane przy użyciu promieni gamma. Sterilizirano z gama sevanjem Стерилизовано гамма-излучением. 감마멸균제품			
<b>Mat PET Ti DMS</b>					
en fr de it es pt ελ cs pl sl ru kr	Dimethyl siloxane, Polyethylene terephthalate, Titanium. Diméthyle siloxane, Polyéthylène téréphthalate, Titane Dimethylsiloxan, Polyethylenterephthalat, Titan. Dimetilsilossano, polietilene tereftalato, titanio. Dimetil siloxano, polietileno tereftalato, titanio. Dimetilsiloxano, tereftalato de polietileno, titânio. Διμεθυλοσιλοξάνιο, τερεφθαλικό πολυαιθυλένιο, τιτάνιο. Dimethylsiloxan, polyetylén teraftalát, titan. Dimetylosiloksan, politereftalan etylenu, tytan. Dimetil siloksan, polietilen tereftalat, titan. Диметилсилоксан, полиэтилентерефталат (ПЭТФ), титан. 디메틸실록산, 폴리에틸렌테레프탈레이트, 티타늄				

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