CYCLOMESH

Sterile Parietal Reinforcement Implants with Drug Delivery System
Mesh to be impregnated with Ropivacaine Kabi 10mg/ml sol inj

en	Instructions for use	Page	2	CE marking date: CE pending
fr	Notice d'instructions	Page	7	Date de marquage CE :En attente de marquage CE
de	Gebrauchsanweisung	Seite	12	Datum der CE-Kennzeichnung: CE pending
it	Istruzioni per l'uso	Pagina	17	Data di marcatura CE: in attesa di CE
es	Instrucciones de uso	Pagina	22	Fecha de marcado CE: CE pendiente
nl	Gebruiksaanwijzing	Pagina	27	Datum CE-markering: CE in behandeling

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CE
0476

NOT290 / 240531
Version du 31/05/2024

This release is the last update of the instructions of use and replace the previous edition
Product under regulation 2017/745 (EU).

Mesh to be impregnated with Ropivacaine Kabi 10mg/ml sol inj
Ropivacaine not provided
Sterile Parietal reinforcement implants with Drug Delivery System
SINGLE USE STERILE PRODUCT
INGUINAL and FEMORAL hernias

MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
ALLEE DES ROSES
59117 WERVICQ SUD
FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

CYCLOMESH

The parietal reinforcement mesh is sterile, non-absorbable and coated with bioresorbable cyclodextrin polymer. This coating has the capacity to adsorb a specific analgesic molecule, Ropivacaine hydrochloride 10mg / ml sol inj (Fresenius Kabi recommended). When the implant is impregnated with Ropivacaine hydrochloride 10mg / ml, the cyclodextrin captures the molecules of Ropivacaine hydrochloride which are then released from the implant.

INDICATIONS

The implant is designed for the treatment of inguinal and femoral hernias. From a clinical point of view, the implant contributes to the treatment of the hernia and the reduction of symptoms such as discomfort and pain. The addition of a pain medicine to the device reduces immediate post-operative pain.

CONTRAINDICATIONS

Do not implant in the following cases :

- Allergy or hypersensitivity to any of the components
- Infected site
- Pregnancy
- Allergy to corn starch and its derivatives.
- Children
- Patient with anticoagulant therapy

Do not impregnate with any molecule other than Ropivacaine hydrochloride 10mg / ml inj sol.

Contraindications and precautions regarding the use of Ropivacaine Kabi 10mg / ml apply and should be taken into consideration when using this device.

TARGET POPULATION

Adult with inguinal or femoral hernia.

IMPLANTED MATERIALS

Implant (non-resorbable): polyethylene terephthalate (multifilament knit).
Coating (resorbable): Cyclodextrin (Poly (citric acid + hydroxypropyl-beta-cyclodextrin))
SILICONE MED-6608-2 (Black silicone ink) according to the references.

Neither human nor animal origin

PERFORMANCES

It has the required mechanical strength to resist to abdominal pressure (16N/cm) and therefore to contribute reinforcing the weak abdominal wall. Once implanted, the macroporous, biocompatible and partially absorbable parietal reinforcement mesh has a mode of action that causes reactive fibrosis leading to its integration and colonization by host tissues

The device also has the capacity to absorb and then locally release Ropivacaine hydrochloride (10mg / ml). The Ropivacaine loading is performed by the surgeon intraoperatively, by impregnating the implant with a Ropivacaine solution (see the procedure for impregnation below). A first analgesic release takes place within the first 5 minutes after implantation (70% of the Ropivacaine loaded), then the rest of the analgesic is released gradually.

The cyclodextrin coating has shown in pre-clinical in vivo tests its ability to absorb Ropivacaine and release it gradually. The claim for reduction of post-operative pain when using the implant with Ropivacaine has been demonstrated by a prospective multicenter randomized clinical study.

LINK TO THE SAFETY AND CLINICAL PERFORMANCE SUMMARY :

https://www.cousin-biotech.com/d-62c425f711511_summary-of-safety-and-clinical-performance-of-sterile-parietal-reinforcement-implants-with-drug-delivery-system.pdf

PRECAUTIONS FOR USE

The implant is not intended to be used only for pain treatment.

The implant should always be separated from the abdominal cavity by the peritoneum.

CYCLOMESH implants are delivered sterile (sterilization with ethylene oxide).

Before any use, check the integrity of the device and the packaging (including the aluminium pouch and the blister). Do not use if the device and / or the packaging is damaged.

Do not use if the device is out of date.

This device must be implanted only by a qualified surgeon trained to the use of the product (knowledge of anatomy and visceral surgery).

The implant is made with polyethylene terephthalate and non metallic component. Consequently, the CYCLOMESH is considered safe for magnetic resonance imaging (MRI safe).

Cousin BIOTECH does not recommend any specific brand of fixing. However, according to the risk analysis, it is recommended to place the fixation points 1cm from the edge of the implant and 1cm apart between each fixation point. Tackers with a square section design are not recommended.

The CYCLOMESH implants are cut-able apart from product references VCBCYFLAPS and VCBCYFLAPM.

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects which include:

- | | | | | |
|--|---------------|------------|--------------------|-------------|
| -Discomfort/pain | -Recurrence | -Adhesions | -Damage to nerves | -Lymphocele |
| -Infection | -Inflammation | -Fistula | -Implant migration | |
| -Allergy to materials and / or coating | -Seroma | -Erosion | -Hematoma | |

The undesirable side effects of Ropivacaine should also be taken into consideration (see the instructions for use for Ropivacaine)

Note: These undesirable side effects may require a second surgery or correction of the original procedure

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss. The surgeon shall also inform the patient that the summary of the safety and performance characteristics of the device and the instructions for use can be found on the COUSIN BIOTECH website.

IMPORTANT: DO NOT REUSE – DO NOT RESTERILIZE

As specified on the labelling, the implant is for single use only. It must not be reused or re-sterilized (potential risks loss of sterility of the implant, risk of infection, loss of efficacy, recurrence).

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.


INSTRUCTIONS FOR IMPREGNATION OF THE MESH

1/ Preparation of the implant with Ropivacaine hydrochloride solution for injection 10 mg / ml (Fresenius-Kabi recommended)

It was designed to be impregnated with Ropivacaine hydrochloride solution for injection 10 mg / ml (Fresenius-Kabi recommended) before implantation:

- Open the box, then the aluminium pouch: it contains a sterile blister
- Place the blister on the sterile operating field
- Follow the instructions printed on the blister's lid:
 - o Inject the required amount of Ropivacaine hydrochloride through the target (see instructions printed on the blister's lid)
 - o Leave to soak for the time indicated on the lid
 - o Open the blister and remove the implant from the blister using atraumatic forceps
 - o Implant the activated mesh


Example of instructions printed on the blister's lid:

Quantity of Cyclodextrin and Ropivacaine per implant:

Lichtenstein							
Shape	Size (cm) Surface (cm ²)	Reference	Quantity of Ropivacaine to inject for impregnation	Time of impregnation	Quantity of Ropivacaine per implant (mg ± 1.1 mg)	Quantity of Cyclodextrin polymer (g ± SD*)	Quantity of PET (knitting) (g ± SD**)
	6x11 57.3	VCBCY0611U	Add 20 mL Ropivacaine chlorhydrate 1% (10 mg/ml)	10 min ≤ ⌚ ≤ 1h	17.8	0.09	0.22
	6x13 67.6	VCBCY0613U			19.3	0.10	0.26
	7x14 87.3	VCBCY0714U			19.4	0.13	0.33
	8x12 84.9	VCBCY0812U			22.6	0.13	0.32
	8x12 99.2	VCBCYFLAPS			22.7	0.15	0.38
	9x14 127.5	VCBCYFLAPM			28.2	0.19	0.48
	14x14 195.1	VCBCY1414U			39.0	0.29	0.74
	14x20 279.1	VCBCY1420U			41.3	0.42	1.06
Laparoscopy							
Shape	Size (cm) Surface (cm ²)	Reference	Quantity of Ropivacaine to inject for impregnation	Time of impregnation	Quantity of Ropivacaine per implant (mg ± 1.1 mg)	Quantity of Cyclodextrin polymer (g ± SD*)	Quantity of PET (knitting) (g ± SD**)
	12x16.5 144.5	VCBCY1216U	Add 20 mL Ropivacaine chlorhydrate 1% (10 mg/ml)	10 min ≤ ⌚ ≤ 1h	31.4	0.22	0.55
↔ Ink marking	SD*: 0.001 g/cm ²						SD**: 0.0005 g/cm ²

Ropivacaine Hydrochloride

See the Ropivacaine instructions for use for the maximum daily dose.

The maximum recommended dose of cyclodextrin is 20mg / kg / day. In the case of the larger implant 14x20cm, containing 420 mg of cyclodextrin polymer, a 70kg patient receives a maximum of 6 mg / kg / day of cyclodextrin polymer, which is less than the maximum recommended amount of 20mg / kg / day.

2/ Implant size

The implant should be of sufficient size with an overlap to cover the hernial orifice beyond the edges of the hernial orifice. Using an implant that is too small can put excessive stress on the fixation points or lead to recurrence of the hernia.

3/ Deployment and fixation

The implant loaded with Ropivacaine hydrochloride 10 mg / ml is placed and then fixed according to the surgeon's usual protocol. Proper fixation is important to secure the placement of the implant.

DESCRIPTION OF THE LICHTENSTEIN TECHNIQUE

- Impregnation of the implant with Ropivacaine: prepare the implant following the operating procedure on the blister's lid (see above)
- The patient is placed in the supine position
- Anesthesia
- Iliac incision of about 5 cm
- Dissection until identification of the hernial sac and identification of nerves
- The hernial sac is identified and reduced
- The transversalis fascia can be tightened using a non resorbable suture. In the event of a direct inguinal hernia, this makes it possible to avoid any recurrence immediately after the operation
- The prosthesis impregnated with Ropivacaine is placed gradually:
 - start by placing the implant on the pubic spine
 - the upper part of the prosthesis is then placed on the joint tendon
 - continue by placing the prosthesis towards the deep inguinal ring, so that the first arm of the prosthesis is placed behind
 - finally place the second arm, so that it surrounds the spermatic cord
 - fix the implant according to the surgeon's habit
 - the apex attaches to the pubic spine
 - the lower edge of the mesh is fixed to the free edge of the inguinal ligament
 - the fixation extends to the anterosuperior iliac spine
 - fix the two ends of the mesh around the spermatic cord
 - the infero-medial corner of the mesh is fixed and overlaps the spine of the pubis
 - fix the mesh to the conjoint tendon
 - the fixation points should be placed 1cm from the edge of the implant and 1cm apart between each fixation point.
 - the cord rests on the prosthesis
 - a suction drain can be placed under the aponeurosis of the external oblique, more particularly in the event of large inguinal hernias, where a significant dissection has been performed
 - the aponeurosis of the external oblique is then closed
 - during perioperative care, prophylactic antibiotics are usually given to the patient within 48 to 72 hours after the procedure

DESCRIPTION OF THE LAPAROSCOPIC TECHNIQUE

1/ Example of repair by TAPP (trans-abdominal-preperitoneal approach technique)

General anaesthesia

Impregnation of the implant with Ropivacaine: prepare the implant by following the operating procedure on the blister's lid (see above)

A pneumoperitoneum is created using the usual technique; intra-abdominal explorations are performed [intra-abdominal pathology, inguinal hernia defects and bags].

The two additional trocars are inserted in direct vision.

The dissection begins and is performed following the classic TAPP technique

The sac of the indirect inguinal hernia should be carefully dissected from the spermatic cord.

The implant is inserted into the intra-abdominal cavity via the trocar, before being deployed over the inguinal region. The marking on the prosthesis is the anatomical landmark for the pubis. The implant is fixed according to the surgeon's usual technique. The operator should avoid the Triangle of Death (or Doom), the Triangle of Pain, the iliac vessels, and should place the fixations laterally to the inguinal ring. The fixation points should be placed 1cm from the edge of the implant and 1cm apart between each fixation point. Tackers with a square section design are not recommended.

The peritoneum should be meticulously closed and no defect should be left open between the peritoneum and the abdominal wall. In addition, it must cover the entire implant.

The additional trocars are removed under direct vision. The fascia at the sub umbilical trocar site should be closed properly.

2/ Example of repair by TEP (extraperitoneal approach technique)

General anaesthesia

Impregnation of the implant with Ropivacaine: prepare the implant by following the operating procedure on the blister's lid (see above)

A 10 mm trocar is inserted. Create the pneumo-pre-peritoneum. Two 5 mm trocars are inserted in direct vision. If a peritoneal defect is created at this time, it should be closed immediately.

The anatomical landmarks are identical to those of a TAPP repair.

The sac of the indirect inguinal hernia should be carefully dissected from the spermatic cord. Direct hernial sacs are easily dissected. It is enough to pull lightly, using an atraumatic forceps. The implant is inserted into the preperitoneal space, before being deployed to the inguinal region. The marking on the prosthesis is the anatomical landmark for the pubis.

The implant is attached to the upper peripheral ligament, and to the upper part of the pubic branch and then to the posterior part of the white line. Fixation continues around and laterally to the inferior epigastric vessels. A few fixations are used to hold the implant laterally to the inner ring. The fixation points should be placed 1cm from the edge of the implant and 1cm apart between each fixation point. Tackers with a square section design are not recommended.

The implant should generously cover the iliac vessels, without straining and without major gaps.

The operator must check that the implant is well anchored to neighbouring structures. The mesh must not migrate and must remain in place. The preperitoneal space should be checked, to avoid possible defects of the peritoneum. If so, they must be closed. The 5 mm trocars are removed in direct vision. The preperitoneal space will collapse. If you suspect that there is a peritoneal defect that

you may have missed, a laparoscopic check-up should be done at the end of the procedure. The trocars are removed, and the fascia defect is closed with the appropriate sutures. The edges of the skin are brought together using the usual technique.

EXPLANTATION AND DISPOSAL OF DEVICES

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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