# LIFT<sup>®</sup>Kit

## STERILE UROGENITAL IMPLANTS: SUBURETHRAL SUPPORT TAPE

en	Instructions for use <b>Surgical Procedure</b>	Page	2 28
fr	Notice d'instructions <b>Technique Opératoire</b>	Page	4 28
de	Gebrauchsanweisung <b>Operationstechnik</b>	Seite	6 28
it	Istruzioni per l'uso <b>Técnica Operatoria</b>	Pagina	8 28
es	Instrucciones de uso <b>Técnica Operatoria</b>	Pagina	10 28
pt	Instruções de uso <b>Técnica Operatoria</b>	Pagina	12 28
ελ	Οδηγίες χρήσης <b>Χειρουργική Διαδικασία</b>	σελίδα	14 28
hu	Használati útmutató <b>Sebészeti eljárás</b>	oldal	16 28
nl	Gebruiksaanwijzing <b>Chirurgische procedure</b>	Page	18 28
pl	Instrukcja obsługi <b>Technika operacyjna</b>	Strona	20 28
sk	Návod na použitie <b>Chirurgický postup</b>	Strana	22 28
cs	Návod k použití <b>Chirurgický postup</b>	Strana	24 28
ro	Instructiuni de utilizare <b>Procedura chirurgicala</b>	Pagina	26 28

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**NOT219/240711**

Version du 11/07/2024

This release is the last update of the instructions for use and replaces the previous edition

## LIFT®Kit Suburethral support tape

### **DESCRIPTION**

The LIFT® kit comprises a tape and a winged guide, 3 ancillary devices (1 pair for the transobturator technique and 1 for the retropubic and suprapubic approach).

They are single use and delivered sterile.

The LIFT® prosthesis consists of a tape of knitted monofilament polypropylene, packaged in a plastic sheath.

It is prolonged by self-tightening loops for placing.

The ancillary devices have different shapes depending on the chosen approach: 1 curved needle for the retropubic and suprapubic approach, a pair of helical needles for the in-out and out-in transobturator approaches.

The end of the ancillary needles has a slot for fixing self-tightening loops.

A single-use, sterile guide is supplied with the implant to secure the passage of the ancillary device during the in-out transobturator approach.

When placing the self-tightening loop in the slot, always make sure it is as tight as possible and that the thread is in the slot designed for this purpose so that the prosthesis does not come apart from the ancillary when passing in the patient's tissues. (See figure 1 below)


Fig 1

### **MATERIALS**

- **Prosthesis:** polypropylene (knitted monofilament)
  - **Thread:** Polyethylene terephthalate + biocompatible dye + dimethyl siloxane
  - **Ancillary device:** Stainless steel needle, polyoxymethylene handle.
  - **Guide:** Liquid crystal polymer + biocompatible dye.
- Not human nor animal origin – Non absorbable.**

### **INDICATIONS**

The LIFT® prosthesis is designed for use as a suburethral support tape for the surgical treatment by vaginal approach in women.

### **PERFORMANCE**

These biocompatible, not biodegradable prostheses cause reactive fibrosis that produces additional urethral support.

These prostheses are easy to insert due to their specially designed shape, excellent flexibility and optimum, rapid colonization and integration.

### **CONTRAINDICATIONS**

Do not implant in the following cases:

- Allergy to any one of the components
- For women receiving anticoagulant therapy.
- Urinary tract infection
- Pregnant women
- Skeletally-immature children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

### **SIDE EFFECTS**

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, urethra, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include hematoma, acute or chronic pain, scarification of the vagina, retraction or release of the mesh, infection, erosion, de novo urgency or voiding difficulties, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

### **USAGE PRECAUTIONS**

The LIFT®Kit device should be implanted by qualified surgeons only (knowledge of gynaecological and urological anatomy and surgery) who are extremely familiar with the device.

The LIFT® Kit implant, guide and ancillary devices are delivered sterile (ethylene oxide sterilization). Check that the packaging is intact before using (do not use if the labels and/or protective packaging are damaged). Do not use if the device is damaged.

The plastic sheath that covers the sling must be removed and must not be implanted.

Do not reuse or resterilize the implant, the guide and ancillary devices.

**IMPORTANT: DO NOT REUSE – DO NOT RESTERILIZE**

As indicated on the product labelling, the implant, guide and ancillary devices are designed for single use. They should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

**STORAGE**

To be stored in a dry place away from sunlight, at room temperature, in its original packaging.

**IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

**EXPLANTATION AND DISPOSAL OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 “Implants for surgery – Retrieval and analysis of surgical implants” Part 1: “Retrieval and Handling”. Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

**INFORMATION REQUESTS AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

**IMPLANT CARD**

**Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.**

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
Installation of the patient, anaesthesia and infiltration performed as specified on the "Preliminary Information" page

1. Bilateral passage of the curved ancillary in the Retzius space to position the LIFT® support tape.
2. Cystoscopic control is done after having passed the cords on each side with out having passed the tape into the tissues.
3. Positioning the tape.
4. Adjustment of the tension can be done by a coughing effort.
5. Cut the tape close to the skin; Homeostasis is checked after adjustment of the tape.

**Suture the vaginal incision with absorbable thread and suture the cutaneous incisions.**

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
Installation of the patient, anaesthesia and infiltration performed as specified on the page named «Preliminary Information 1-2-3»

1. Trans-obturator incisions of approximately 1 cm each, positioned at least 1 cm above the horizontal line passing by the urethral meatus at 1 cm from the ischiopubic branch.
2. Vaginal incision and urethro-vaginal dissection.
3. Passage of helicoid ancillary by way of the trans-obturator (blue handle). (the tunnelling device is inserted into one of the trans-obturator incisions).
4. • Hooking of LIFT® bandage onto needle.
5. • Passage of the tape; and heterolateral passage: performance of the 3 previous operations on the heterolateral side.
6. • The tape is deployed by pulling on the thread. There should remain a space of about one centimetre between the tape and the urethra.
7. • Cutting of the tape, flush with the skin at the level of the incisions. Hemostasis is checked after tensioning of the tape.

• Closure of the vaginal wall by a crossed suture with absorbable thread and closure of the cutaneous incision.

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