NAJATM

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NOT238 / 240711

Version du 11/07/2024

This release is the last update of the instructions of use and replace the previous edition

NAJA ancillaries set

Ligament system for spinal stabilisation

DESCRIPTION

This ancillaries set is intended for the placement of the NAJA ligament system for spinal stabilisation. The ancillaries are Class I medical devices, supplied non-sterile, intended for temporary use and reusable.

MATERIALS

The different elements of the set are made from:

- Titanium
- Stainless steel
- Acetal copolymer
- Radel®
- Medthin 43
- Silicone

All the materials used are non-absorbable and are neither of human nor animal origin.

PRECAUTIONS FOR USE

The ancillaries must be stored in a container intended for that purpose or in equivalent packaging. Do not let the instruments come into contact with objects which could alter the surface. Carry out a visual examination of each instrument before each use in order to detect any possible deterioration. The ancillaries must not be used except with elements corresponding to the original COUSIN BIOTECH system. Use of the specific ancillary COUSIN BIOTECH system is required for implantation of the implant. In addition, the instruments must not be used except by a surgeon with appropriate training in the placement of NAJA implants.

DECONTAMINATION AND CLEANING OF ANCILLARIES

Decontamination and cleaning are performed under the responsibility of the healthcare centers. The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011). This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments

The detachable instruments have to be dismantled

☐☐ Products advised for the cleaning

- an enzymatic detergent neutral
 - Or products adapted for the cleaning (neutral or soft alkaline) of chirurgical instruments in compliance with the regulations in effect.
- **Warning:** Don't use corrosive factor or caustic cleaning product

2 process are possible

A- Chart Automated cleaning in a disinfectant washer with manual Pre-cleaning:

| | Step | Step Description | Step Instruction | Accessories | Duration |
|----------------|------|------------------------|--|--|---|
| Cleaning Steps | 1 | Contaminati on Removal | Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush | - Tap water - Soft bristle brush (Do not utilize metal cleaning brushes) | Until all visible soi lis removed |
| | 2 | Pre-soak | Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position. | - Detergent - Water jet pistol | Until product is visually clean |
| | 3 | Ultrasonic cleaning | Clean the instrument in the ultrasonic tray including the solution. Be careful every unclench connector or screws and bolts by the vibrations will be reassure. | - Ultrasonic tray | Minimum duration 10 minutes, this duration can vary according product |
| | 4 | Soak | Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position. | - Reserve osmosis water | Minimum 30 seconds of soak |
| | 5 | Automated Washer | Place entire disassembled device into the automated washer | - Automated Washer - Disassembled Device | Minimum total cycle time: 39 minutes |

| Instructions of cleaning in washer / disinfectior | | | |
|---|------------|---------------------|-------------------------|
| Cycles parameters | Time | Minimum temperature | Type of detergent/water |
| Pre-cleaning | 2 minutes | <45°C / <113°F | Tap water |
| Cleaning | 5 minutes | 55°C / 131°F | Detergent |
| Draining | 2 minutes | <45°C / <113°F | Tap water |
| Thermic draining | 5 minutes | 90°C / 194°F | Reserve osmosis water |
| Drying | 25 minutes | > 70°C / > 158°F | Not applicable |

| | | | | | |
|-----------------|---|--------------|--|--------------------------------|-------------------------------|
| Finishing Steps | 6 | Final Rinse | Rinse carefully with distilled water | - Distilled water | Minimum 1 minute |
| | 7 | Final Drying | Dry devise utilizing medical quality filtered air | - Medical quality filtered air | Until product is visually dry |
| | 8 | Tidying | The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization. | | |

B- Chart Just manual cleaning :

| | Step | Step Description | Step Instruction | Accessories | Duration |
|--------------------|------|--------------------------|--|--|--|
| Cleaning Steps | 1 | Contamination Removal | Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush | <ul style="list-style-type: none"> - Tap water - Soft bristle brush (Do not utilize metal cleaning brushes) | Until all visible soil is removed |
| | 2 | Drying | Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available | <ul style="list-style-type: none"> - Non-scheduling wipe - Medical quality filtered compressed air | Until product is visually dry |
| Disinfection Steps | 3 | Disinfection Application | Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument. | <ul style="list-style-type: none"> - Detergent | Minimum duration 1 minute, this duration can vary according product |
| | 4 | Drying | After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water. | <ul style="list-style-type: none"> - Disposable cloth - Detergent | Until product is visually dry |
| | 5 | Manual Disinfection | Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution. | <ul style="list-style-type: none"> - Soft bristle brush - Detergent - Syringe | Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned |
| Disinfection Steps | 6 | Rinse and finish | Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water | <ul style="list-style-type: none"> - Reverse osmosis water - Syringe | Minimum 1 minute |
| | 7 | Final Rinse | Rinse carefully with distilled water | <ul style="list-style-type: none"> - Distilled water | Minimum 1 minute |
| | 8 | Final Drying | Dry device utilizing medical quality filtered air | <ul style="list-style-type: none"> - Medical quality filtered air | Until product is visually dry |
| | 9 | Tidying | The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization. | | |

Note: In case of **suspicion** of **Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

STERILIZATION OF ANCILLAIRES

Sterilization is performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

1- The sterilization in autoclave is to make in a **specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

NB: It **isn't recommended** to realize the sterilization with peelables pouches **autoclavable** without a container □ Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

| Steam sterilization parameters | Values | Values |
|---|--------------------|--------------------|
| Cycle Type | Pre-vacuum | Pre-vacuum |
| Set Point Temperature | 134°C / 273°F | 132°C / 269.6°F |
| Cycle Time | Minimum 3 minutes | Minimum 4 minutes |
| Dry Time | Minimum 20 minutes | Minimum 20 minutes |
| Cool Time (inside and outside the autoclave chambers) | Minimum 40 minutes | Minimum 40 minutes |

3 - Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrity and the correct operation, of the locking systems (clipping, prehension) without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent sterilization process for which it will be fully responsible.

REQUESTS FOR INFORMATION AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to doing everything possible to produce and supply a quality medical device. However, in the event that a health professional (client, user, prescriber, etc.) has a complaint or a reason for dissatisfaction regarding a product, in terms of quality, safety or performance, they must inform COUSIN BIOTECH of it as soon as possible.

In case of an instrument malfunctioning or if an instrument has contributed to causing a serious undesirable effect for the patient, the health centre must follow the legal procedures in force in its country and must inform COUSIN BIOTECH of this as soon as possible.

For all correspondence, please specify the reference, the batch number and the contact details of a referent, as well as a comprehensive description of the incident or complaint.

The brochures, documentation and operating technique are available through a simple request to COUSIN BIOTECH and its distributors.

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