BIORING®

ADJUSTABLE GASTRIC BANDING & PORT

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	5
de	<u>Gebrauchsanweisung</u>	Seite	8
it	<u>Istruzioni per l'uso</u>	Pagina	11
es	<u>Instrucciones de uso</u>	Pagina	14
ελ	<u>Οδηγίες χρήσης</u>	Σελίδα	17
ro	<u>Instrucțiuni de Utilizare</u>	Pagină	21
cs	<u>Pokyny k použití</u>	Strana	24
et	<u>Kasutusjuhend</u>	lehel	27
lt	<u>Naudojimo instrukcija</u>	Puslapis	30
lv	<u>Lietošanas instrukcija</u>	Lappuse	33
nl	<u>Gebruiksaanwijzing</u>	Pagina	36

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Gastric ring kit + port

DESCRIPTION

The **BIORING®** adjustable gastric band is composed of:

- A radiopaque body, comprising an external reinforcement and an inflatable inner part. The extremities of the body form the closure system of the band after its placing around the stomach.
- A radiopaque injection port, which is usually placed in the left subcostal or presternal region of the patient. This site enables the percutaneous adjustment of the diameter of the adjustable gastric band.
- A radiopaque catheter, which links the body of the adjustable gastric band to the injection port. The ensemble is non-absorbable.
- A stiffener, which secures the connection between the injection port and the catheter via a click-shut system.

COMPOSITION

The body is composed of dimethylsiloxane. The external reinforcement is made of polyester which is fully impregnated with dimethylsiloxane, barium sulphate, and a biocompatible dye.

The catheter is made of dimethylsiloxane impregnated with barium sulphate, and printed with biocompatible ink.

The injection port is composed of polyoxymethylene impregnated with barium sulphate. Its septum is composed of dimethylsiloxane.

The stiffener is made of dimethylsiloxane and barium sulphate.

No human or animal origin. Non-resorbable.

INDICATIONS

The **BIORING®** adjustable gastric band system is indicated for the reduction of excess body weight in persons affected with morbid obesity, and who are resistant to dieting, i.e. persons whose Body Mass Index (BMI, formula = weight in kilograms / (height in metres)²) is greater than or equal to 40 kg/m², or 35 to 40 kg/m² when combined with at least one comorbidity. The surgeon is free to prescribe the adjustable gastric band at his own discretion.

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with dimethyl siloxane and non metallic component which are non-ferromagnetic in these quantities according to literature. Moreover, the implant is presented as MRI safe in the misafety.com database with 1,5 and 3 Tesla. Consequently, the Bioring® is considered MR safe.

CONTRAINDICATIONS

- Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract.
- Severe organic disease.
- Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs.
- Diagnosis or pre-existing symptoms or family history of connective tissue autoimmune disease.
- Pregnancy.
- Patients with a drug or alcohol addiction.
- Patients who refuse to accept the dietary restrictions required by the procedure.
- Serious psychiatric disorders.
- Do not implant the device at the site of infection or in a septic environment.
- Do not implant the device in the event of allergy to the components.
- Do not implant the device in growing children.

POSSIBLE COMPLICATIONS

Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, the surgical technique, and the degree of intolerance of the patient to any foreign body implanted in the body.

Perioperative complications (due to the surgical technique):

- Lesions of the spleen
- Perforation of the stomach or oesophagus
- Other bleeding

In the immediate postoperative period:

- Nausea, vomiting (especially in the first few days after the intervention, and when the patient eats more than recommended).
- Bloating due to gas
- Dysphagia, dehydration
- Thromboembolic complications
- Infectious complications
- Postoperative bleeding

Delayed postoperative complications

- Band dilation
- Band erosion
- Oesophageal dilation
- Pouch dilation
- Band alteration
- Band intolerance
- Pain located at the site of port implantation
- Band infection
- Slipping of the band,
- Gastroesophageal reflux,
- Pyrosis,
- Gastritis of the upper pouch,
- Migration of the band,
- Displacement of the injection port,
- Infection of the injection port,
- Constipation

- Constrictive obstruction (due for example to oedema, food, an inadequate calibration, slippage of the band, or to torsion or dilation of the upper pouch),
- Deflation of the band (due to a leak in the band, injection port, or connection tube),
- Rapid weight loss may lead to symptoms of malnutrition, anaemia, and associated complications

WARNING

The laparoscopic implantation technique is a state-of-the-art procedure. Surgeons who intend to perform this procedure should possess the necessary equipment and have extensive experience in laparoscopy.

PRECAUTIONS

- Store in a dry place, away from sunlight and at room temperature in its original packaging.
- Check the integrity of the packaging (do not use if the labels and / or protective packaging are damaged).
- Do not use if the device is damaged or has past its expiry date.

DO NOT REUSE

DO NOT RE-STERILIZE

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be performed in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis, following the current protocol. This protocol is available on demand from COUSIN BIOTECH. It is important to note that any implant that is not to be cleaned and disinfected before expedition must be contained in a sealed package. Explanted medical devices must be eliminated in accordance with the standards applicable in the country for the disposal of infectious hazardous waste. There are no specific recommendations for the elimination of a non-implanted device. The second paragraph is in case of an information request or complaint.

INFORMATION REQUEST AND COMPLAINT

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, he/she must inform COUSIN BIOTECH of this as soon as possible. In the case of failure of an implant or if it contributed to causing serious adverse effects for the patient, the health centre must follow the legal procedures in its country and inform COUSIN BIOTECH of this as quickly as possible. For any correspondence, please quote the reference and batch number, the reference coordinates and give a comprehensive description of the incident or complaint. Brochures, documentation and the surgical technique are available on request from COUSIN BIOTECH and its distributors.

SURGICAL TECHNIQUE

Material provided

- | | |
|-----------------------------------------------------------------------------------|-------------------------------------------|
| - One adjustable, radiopaque gastric band, with closure system | - One radiopaque injection port |
| - One syringe | - One radiopaque catheter flushing device |
| - One Huber needle | - One handling tag |
| - A gastric tube is available separately for the calibration of the small stomach | - One radiopaque stiffener |

Preparation of the material: *All of the procedures should be performed in a sterile field.*

- Flush the band with physiological saline using the syringe and catheter rinsing device provided to completely evacuate the air from the device.
- Inspect the inflatable part of the band to detect any irregular inflation or leaks. When handling the band, make sure that it does not come into contact with sharp or pointed instruments.
- Clamp the catheter to prevent air from entering the device again.
- Place the handling tag at the free end of the catheter.
- Flush the implantable injection port with physiological saline, using the syringe and Huber needle.

Furthermore, to get a maximal internal diameter of the band, it is advisable to aspirate the content of the adjustable gastric band completely, using the syringe and catheter-rinsing device, especially during procedures using the "pars flaccida" technique, and in patients with a significant quantity of perigastric adipose tissue.

Surgical technique:


Prophylactic Antibiotics: It is recommended to deliver an antibiotic prophylaxis against skin and gut flora, prior to beginning the procedure.

Various surgical techniques may be used, including: the perigastric technique, the "pars flaccida" technique, or a combination of these two techniques.

Adjustable gastric bands are usually placed laparoscopically, under general anesthesia:

- The anaesthetist introduces the cuffed gastric tube into the stomach, making it possible to calibrate the small stomach.
- Dissect the lesser curvature or the right crus.
- 2nd window next to the His angle.
- Retro-gastric tunnelling.
- Introduction of the adjustable gastric band into the abdomen. When passed through the trocar (12mm minimum), the clip must be positioned on the reinforced part of the buckle.
- Inflate the cuff of the gastric tube with 15 cm³ (ml) of physiological saline.
- Correctly position the band.

- Suture the band, using 3 to 6 sutures, to the gastric serosa, upstream and downstream (anterior gastro gastric suture enclosing the band in a "tunnel").
- Connection of the tube to the injection port.
- General port placement recommendations:
 - Make sure the catheter shows no sharp bending
 - When connecting the catheter to the port make sure you don't create twisting/spring effect on the catheter: check that the port remains totally flat on a surface
 - Make sure the stiffener is properly clipped to the port (click-slot system).
 - In the abdominal wall layer, the path of the catheter must be oblique rather than perpendicular to the port. This prevents lifting effect on the port which with time could lead to port rotation:

The injection port may be:

- inserted and sutured in the aponeurosis of the rectus abdominis muscle under the left costal edge,
- inserted and sutured subcutaneously under the left costal edge,
- inserted and sutured pre-sternally, on the xiphoid apophysis;
- inserted and sutured in any other position according to the surgeon's appreciation.

The bands may also be implanted via laparotomy.

The band can be inflated immediately using physiological saline (NaCl 0.9%). ONLY use a **Huber needle** to fill-in or deflate the band. The volume, as well as the interval before filling, is at the discretion of the practitioner/surgeon (maximum volume, see the next table).

Inflating of the band: maximum volume of physiological saline according to sizes of the bands

Band sizes (diameter : mm)	Maximum volume of physiological saline (cm ³)
Ø 21/24	10
Ø 27/29	10
Ø 36/38	15

An example of band-filling algorithm is available on request.

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