# EXTENDED PLATE FOR CORTICAL FIXATION SYSTEMS

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**NOT252\_240919**

Version du 19/09/2024

**Rx Only** Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician

This release is the last update of the instructions of use and replace the previous edition

**EXTENDED PLATE FOR CORTICAL FIXATION SYSTEMS  
STERILE CORTICAL FIXATION DEVICES  
SINGLE USE STERILE PRODUCT**

**DESCRIPTION**

Extended Plate to be used with cortical fixation systems (adjustable and continuous loop) to allow a larger surface of contact with the cortical bone.

**MATERIALS**

**- *Materials implanted:***

Extended plate: Non-ferromagnetic titanium alloy (Ti TA6V ELI)

Neither animal nor human origin. Non-absorbable

**INDICATIONS**

The Cortical Fixation Systems are used for fixation of tendons and ligaments during orthopaedic reconstruction procedures for the anterior cruciate ligament (ACL) reconstruction.

**PERFORMANCE**

This system facilitates ligament or tendon repair/reconstruction by graft fixation in the bone.

It has been approved for long-term implantation in the human body.

**CONTRAINDICATIONS**

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Insufficient quantity or quality of bone.
- Blood supply and previous infections which may tend to retard healing.
- Active infection.
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- Not to be used for the bone – patella tendon – bone technique.
- Do not implant in children while they are growing.

**POSSIBLE ADVERSE EFFECTS**

Possible complications linked to the implant:

- 1- Breaking of the loop, breaking of the traction wires, breaking of the button
- 2- Failure when placed under tension
- 3- Migration or upturn of the button
- 4- Reaction to a foreign body

Possible complications linked to insertion of the implant:

- 1- Pain in the incision area or surgical site
- 2- Local and transient inflammatory reactions
- 3- Infection
- 4- Wound dehiscence

**WARNING**

- 1- Any decision to remove an implant must take into account the fact that a second surgical operation may involve a risk for the patient. Removal must be followed by suitable postoperative care.
- 2- The patient must be informed of the limitations inherent in the nature and function of the implant, in particular regarding the weight and constraints placed by the body on the system before the bone has sufficiently healed.
- 3- Ensure that the bone tunnel has the correct length and diameter.

**PRECAUTIONS BEFORE USE**

Before use, check the integrity of the implant and the packaging.

Do not use if the implant, labels and/or protective packaging is damaged.

Do not use the product after the expiry date indicated on the label. This could compromise the performance of the system and the patient's safety.

Unpacking and handling of the implant must be carried out in accordance with asepsis standards, to ensure sterility.

When opening the cardboard box, the surgeon must check that the batch number and model match to those indicated on the label on the inside cover. If they do not, use of the implant is strictly forbidden.

The device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

This implant can be used with standard instruments available on the market and specified for orthopaedic surgery for ligament or tendon anchorage.

Surgical instruments are subject to wear during normal use. Instruments that have been used many times or been subject to very high stress are likely to break. Instruments must only be used for their intended purpose. COUSIN SURGERY recommends that instruments should be regularly inspected for wear and bending.

**DO NOT REUSE**

**DO NOT RESTERILIZE**

The implant is for single use, as specified on the label. It must not on any account be reused or re-sterilized (potential risks include, but are not limited to: loss of sterility of the product, risk of infection, loss of effectiveness of the product).

COUSIN SURGERY does not provide any guarantee or recommendations regarding the use of a specific make of fixation device. The properties of these devices are subject to alteration by the manufacturer, over which COUSIN SURGERY has no control.

This device must **only** be implanted by a qualified surgeon trained in the use of the product. Avoid using staples or sharp-edged screws.

## **SURGICAL TECHNIQUE**

### **Preliminary information**

The system is adapted for the use of the right internal and semitendinosus tendons with the all-inside and inside-out techniques as well as retrodrilling. **Adjustable cortical fixation systems:** In case of retrodrilling, COUSIN SURGERY recommends a tunnel diameter of 4.5 mm to insert the implant. In case of a full tunnel, it is preferable to use the additional extended plate.

Use of the device with a tunnel diameter of less than 4.5 mm has not been validated. If the tunnel diameter is greater than 4.5 mm, the surgeon can decide whether or not to use the extended plate.

**Cortical Fixation Systems:** In case of retrodrilling, COUSIN SURGERY recommends a tunnel diameter of 5 mm to insert the implant. In case of a full tunnel, it is preferable to use the additional extended plate.

Use of the device with a tunnel diameter of less than 5 mm has not been validated. If the tunnel diameter is greater than 5 mm, the surgeon can decide whether or not to use the extended plate.

### **Technique**

1. Suture the graft and put it in the loop.

#### **2. Adjustable cortical fixation systems:**

- Create a bone tunnel to insert the implant with a diameter of 4.5 mm (recommended diameter) and a graft cavity suited to its diameter.

#### **Cortical Fixation Systems:**

- Create a bone tunnel to insert the implant with a diameter of 5 mm (recommended diameter) and a graft cavity suited to its diameter.

#### **3. Adjustable cortical fixation systems:**

- go directly to step 4.

#### **Cortical Fixation Systems:**

Choosing the size of the implant. Total length of the femoral tunnel (between the two cortices) - length of the graft tunnel +5 cm. If the result is not equal to the size of the implant, round up to the larger implant size.

#### **4. Adjustable cortical fixation systems:**

Insert the suture threads and adjustment threads through the bone tunnel.

#### **Cortical Fixation Systems:**

Pass the threads through the bone tunnel.

5. Pull on one of the sutures until the button comes out of the bone tunnel through the cortical bone.

6. Tip over the button against the cortical bone surface by pulling on the graft.

7. Remove one of the two sutures

8. Slide the extended plate under the button. Ensure that the remaining suture is on the split side to avoid hindering the interlocking.

#### **9. Cortical Fixation Systems:**

- go to step 12 and then step 14.

#### **10. Adjustable cortical fixation systems:**

While still applying counter tension to the graft, slide the adjustment threads one after another to adjust the loop to the required length, starting with the white thread. Do not apply tension to the black thread.

#### **11. Adjustable cortical fixation systems:**

To lock the loop, continue to apply tension to the white adjustment thread and pull firmly on the black adjustment thread.

12. Fix the graft in the opposite bone tunnel using the required method.

#### **Adjustable cortical fixation systems:**

- go to step 13.

#### **Cortical Fixation Systems:**

- go to step 14.

13. If necessary, slide the adjustment threads alternately once again to obtain the required tension.

And repeat step 11.

14. Remove the suture thread(s) from the button.

15. **Adjustable cortical fixation systems:** Section the adjustment threads at a minimum distance of 10 mm from the cortical button.

## **STORAGE PRECAUTIONS**

Store in a dry place away from sunlight, at room temperature and in the original packaging.

## **EXPLANTATION AND DISPOSAL OF THE DEVICE**

COUSIN SURGERY recommends that explantation and handling should be carried out in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

All explanted devices must be returned for analysis according to the protocol in force. This is available on request from COUSIN SURGERY. It should be noted that any implants that cannot be cleaned and disinfected before shipping must be placed in a sealed package.

Disposal of explanted medical devices must be carried out in accordance with the standards in force in the country, concerning the disposal of infectious waste.

There are no special recommendations for the disposal of non-implanted devices.

### **INFORMATION REQUEST AND COMPLAINTS**

In accordance with its quality policy, COUSIN SURGERY undertakes to make every effort to produce and supply a high quality medical device. However, if a health professional (customer, user, consultant etc.) has a complaint or grounds for dissatisfaction with a product, regarding its quality, safety or performance, he/she must inform COUSIN SURGERY as soon as possible.

In the event of malfunctioning of an implant or if it contributes towards a serious adverse effect for the patient, the healthcare centre must follow the legal procedures in force in the country, and immediately inform COUSIN SURGERY.

In all correspondence, please state the reference, batch number, name and address of a consultant, and provide a full description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN SURGERY and its distributors.

### **IMPLANT CARD**

**Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.**

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