SACROMESH® SOFT PROLAPS

STERILE UROGENITAL IMPLANTS FOR THE TREATMENT OF PROLAPSE

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Made in France

NOT307_240215
Version du 15/02/2024

This release is the last update of the instructions of use and replace the previous edition

Product under regulation 2017/745 (EU)

SACROMESH® SOFT PROLAPS

STERILE UROGENITAL IMPLANTS FOR THE TREATMENT OF PROLAPSE
SINGLE USE STERILE PRODUCT

MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
ALLEE DES ROSES
59117 WERVICQ SUD
FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

Sacromesh® Soft Prolaps sterile implants consist of knitted macroporous polypropylene monofilament:

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They are intended for sacrocolpopexy of the pelvic prolapse in women, by abdominal approach.

INDICATIONS

Sacromesh® Soft Prolaps prosthesis are designed for use as reinforcement and support of tissues when treating genital prolapse surgically by sacrocolpopexy.

They can be used in either conventional or laparoscopic surgery.

CONTRAINDICATIONS

Do not use the implant in the following circumstances:

- Allergy to any of the components
- In a place with progressive or latent infection or signs of tissue necrosis
- For women receiving anticoagulant therapy
- Urinary or genital tract infection or non-aseptic conditions
- Pregnancy
- Growing children

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:

- Discomfort / pain (Acute or chronic pain)
- Infection
- Inflammation
- Local irritation
- Foreign body reaction
- Recurrence
- Hematoma
- Fistula formation
- Dyspareunia
- Allergic reaction
- Mesh migration / displacement
- Retraction / Release of the mesh
- Adhesion formation
- Occlusion / Obstruction
- Retention of urine
- Perforation or damage to blood vessels, bladder, urethra, rectum, intraperitoneal organs, or nerves
- Vaginal extrusion
- Erosion (through the urethra or surrounding tissues)
- Scarification of the vagina
- Sexual dysfunction
- Adverse obstetric outcomes (preterm premature rupture of membranes)
- De novo urgency or voiding difficulties

One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

TARGET POPULATION

Adult (grown up) female for the surgical treatment.

IMPLANTED MATERIALS

Polypropylene (non resorbable)

References	Weight of Polypropylene (g)
SACROMESH9	0,56 g +/- 0.1 g
GCBPRANT1U	0,27 g +/- 0.1 g

Origin neither human nor animal. Non resorbable.

PERFORMANCES

The device Sacromesh® Soft Prolaps contributes to treat the genital prolapse and reduces its recurrence. This range is biocompatible designed for a reinforcement and support of tissues with a separate tension adjustment.

LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:

https://www.cousin-biotech.com/d-6318526d8b749_under-construction.pdf

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with polypropylene and non-metallic component. Consequently, Sacromesh® Soft Prolaps is considered MR safe.

PRECAUTIONS FOR USE

Sacromesh® Soft Prolaps meshes are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which peelable pouches). Do not use if the device is out of date.

Do not use in the event of deterioration of the device and/or the packaging.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the mesh are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Prevent the mesh from any contact with objects that may alter their surface.

Sacromesh® Soft Prolaps devices must always be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments, and the surgical technique.

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

IMPORTANT: DO NOT REUSE – DO NOT RESTERILIZE

As specified on the product's labelling, the Sacromesh® Soft Prolaps mesh is for single use only. It cannot be re-used and/or resterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence).

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal.

An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer. Please refer to the section "Implant card information leaflet" for instructions on how to complete it. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.**

The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

STORAGE PRECAUTIONS

Store in a dry place away from sunlight and at room temperature in its original packaging.

SURGICAL TECHNIQUE

- Surgery under general anaesthetic.
- Abdominal incision or installation of trocars.
- Exposure of the lumbosacral angle.
- Dissection of the prerectal fascia after incision of the anterior wall of Douglas' cul-de-sac. The dissection continues laterally down to the posterior elevator muscles, using a vaginal valve if necessary.
- The mesh is fixed as distally as possible on the pelvic floor.
- Dissection of the vagina's anterior wall.
- The anterior mesh is fixed by a running stitch, avoiding transfixation. Depending on the mesh used by the surgeon and of his preference, the mesh is passed laterally to the uterine isthmus. The uterus is held in place, two meshes are passed on either side of it (for precut mesh or mesh cut longitudinally by the surgeon).
- Fastening the anterior and posterior mesh to the Promontory (presacral ligament) without excessive tension.
- Current recommendation is to do this stitch on the Promontory with a non-resorbable suture and avoid the use of staplers or tackers.
- Peritonization must be as complete as possible to avoid contact between the prosthesis and the intraperitoneal constituents.

EXPLANTATION AND DISPOSAL OF DEVICES

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

Carte Implant - Implant Card		Carte Implant - Implant Card	
	A conserver / Keep it		
Patient name		_____	 COUSIN SURGERY <i>We care for Surgery</i>
Date of surgical procedure		. . / . . /	
Name of practitioner		_____ _____	

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
<ul style="list-style-type: none"> • Indication that the product is a medical device • Description of the device family • Trade name of the device • Reference of the device • Device is MR safe • Denomination name of the device in sales countries 	<ul style="list-style-type: none"> • UDI number • Batch number • Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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