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Version du 30/10/2024

www.cousin-surgery.com

Made in France
Caution: Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

Sterile spinal dynamic posterior stabilization device

1-DESCRIPTION

The sterile spinal dynamic posterior stabilization device BDyn is intended to restore the stabilization of the non-cervical segment of the spine by preserving the anatomical lordosis and the deadening of the intervertebral joint. It is composed of the BDyn spinal shock absorber fixed on the vertebrae thanks to the dedicated screws and can be replaced according to the decision of the surgeon by a dedicated replacement rigid rod.

2-MATERIALS

The BDyn Device components are manufactured from medical grade titanium TA6V (ISO 5832-3; ASTM F 136), silicone, long-term implantable polyurethane.

BDyn implant: Titanium alloy Ti6Al4V ELI (ISO 5832-3; ASTM F 136), Polydimethyl siloxane (PDMS) MED 4770 unrestricted Polycarbonate urethane (PCU) Bionate® II 80A

BDyn ancillaries: Stainless steel, Titanium, Polypropylene
No human or animal origin – Non resorbable.

3-INDICATIONS

The sterile spinal dynamic posterior stabilization device BDyn is intended for posterior stabilization from thoracic vertebrae T10 to sacrum S1 with or without bone graft for the following indications:

- Degenerative intervertebral disc disease and/or articular facets confirmed by further examinations
- Spinal canal stenosis
- Degenerative spondylolisthesis grade 1
- Segmental instability

4-PERFORMANCE

The BDyn spinal shock absorber of the sterile spinal dynamic posterior stabilization device BDyn consists of a metallic cylindrical hollow part containing elastomer components made of silicone elastomers (polydimethyl siloxane (PDMS)) and long-term implantable polycarbonate urethane (PCU) which are bending out under the effect of a metallic piston rod connected with the vertebra of the treated segment by the pedicular screws tested and/or approved by the company COUSIN BIOTECH. The elastomer components assure the absorption of the mechanical loads applied on the intervertebral joint in compression, pulling, flexion-traction and lateral flexion. Some configurations enable fixation across several spinal segments.

5-CONTRAINDICATIONS

- Active infectious process or significant risk of infection (immunocompromise)
- Signs of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Mental illness
- Grossly distorted anatomy caused by congenital abnormalities
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of the sedimentation rate unexplained by other diseases, elevation of the white blood count
- Suspected or documented metal allergy or intolerance
- Any case where the implant components selected for use would be too large or too small to achieve a successful result
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance
- Any patient unwilling to follow postoperative instructions
- Any case not describe in the indications
- Traumas (i.e. fracture or dislocation)
- Abnormal curvatures (i.e. scoliosis and/or hyper lordosis)
- Tumors.
- Spondylolisthesis grade 2 and more
- Pseudarthrosis and/or failed previous fusion
- Severe bone resorption, osteomalacia, severe osteoporosis

6-UNDESIRABLE SIDE EFFECTS

All of the possible adverse events associated with spinal surgery and without instrumentation are possible:

- Infection
- Pseudomeningocele, fistula, breach dura, persistent CSF leakage, meningitis
- Loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits
- *Cauda equina* syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss
- Urinary retention or loss of bladder control or other types of urological system compromise
- Scar formation possibly causing by a neurological compromise or compression around nerves and/or pain
- Fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery
- Herniated *nucleus pulposus*, disc disruption or degeneration at, above, or below the level of surgery, canal adjacent stenosis
- Non-union or pseudarthrosis, delayed union. Mal union
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function

- Inability to perform the activities of daily living
- Bone loss or decrease in bone density
- Graft donor site complications including pain, fracture, or wound healing problems
- Ileus, gastritis, bowel obstruction or loss of bowel control or other types of gastrointestinal system compromise
- Hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise
- Reproduction system compromise, sterility, sexual dysfunction
- Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc
- Change in mental status
- Death

All of the possible adverse events associated with spinal surgery with instrumentation are possible. A listing of potential adverse events linked to the medical device includes, not limited to:

- Early or late loosening of any or all of the components
- Disassembly, bending and/or breakage of any or all of the components (screw breakage)
- Foreign body (allergic) reaction to implants, debris, corrosion products (from crevice, fretting, and/or general corrosion), including metallosis, tumor formation and/or autoimmune disease
- Pressure on the skin from components parts with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, fibrosis, neurosis, and/or pain
- Tissue or nerve damage caused by improper positioning and placement of implants or instruments
- Post-operative change in spinal curvature, loss of correction, height, and/or reduction

7-PRECAUTIONS FOR USE

The BDyn devices are delivered sterile, sterilized by ethylene oxide.

Before any use, inspect the integrity of the packaging and device (including peelable pouches).

Do not use in the event of deterioration of the labels and/or the device and/or the packaging.

Do not use if the device is out of date.

The sterile spinal dynamic posterior stabilization device BDyn is only intended to adult human female or male more than 40 Kg.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

The installation of BDyn device can only be made with pedicular screws tested and/or approved by the company COUSIN BIOTECH.

The eventual substitution of the BDyn spinal shock absorber has to be carried out only with the dedicated replacement rigid rods.

You should never use implants made of stainless steel and titanium alloy in the same construct.

A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise the results. This device system is not intended to be the sole means of spinal support. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly and/or breakage of the device(s) will eventually occur.

Preoperative and operating procedures, including precise knowledge of suitable surgical techniques, and proper selection of the good reference of the device adapted to the patient and its narrow setting up are important considerations in the successful use of the device by the surgeon.

The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique

Postoperative precautions

MRI safety

The BDyn implant has been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). Non-clinical testing has demonstrated that the BDyn orthopaedic implants are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions. Failure to follow these conditions may result in injury.

Maximum Static Magnetic Field Strength	1.5T / 3.0T
Type of Nuclei	Hydrogen
Static Magnetic Field (B0) Orientation	Horizontal
Magnet Type	Cylindrical-bore
Maximum Spatial Field Gradient	87.1 T/m (8710 gauss/cm)
Maximum B0*ldB0/drl product	101.5 T2/m
RF Excitation	1.5T: Circularly Polarized (CP) 3T: Multichannel-2 (MC-2)
RF Transmit / Receive Coil Type	Integrated Whole Body Transmit and Receive Coil
RF Power	Normal Operating Mode
Maximum Whole-Body SAR	Whole body SAR ≤ 2 W/kg Note: Under the scan conditions defined above, the BDyn implants are expected to produce a maximum temperature rise of less than 6.1 °C at 1.5T and less than 3.4 °C at 3T after 15 minutes of continuous scanning.
Patient Characteristics	1.5T: Patients with uncompromised thermoregulation and under controlled conditions (a medical doctor or a dedicated trained person can respond instantly to heat induced physiological stress)

	3T: Patients with uncompromised thermoregulation and under uncontrolled conditions, or patients with compromised thermoregulation (all persons with impaired systemic or reduced local thermoregulation) under controlled conditions.
MR Image Artifact	MR image quality may be compromised if the imaging area of interest is in the same area of the implant. Some manipulation of scan parameters may be required to compensate for the artifact. In non-clinical testing, the image artifact caused by the device extends approximately 36 mm from the implant when imaged with a spin echo pulse sequence and 22 mm with a gradient echo, both at 1.5T.

Consequently, the BDyn is considered MR conditional.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

In accordance with the labelling of this product, BDyn implants are for single use; they must not in any case be reused or re-sterilised (potential risks include, but are not limited to, loss of sterility of the products, risk of infection, loss of product efficacy, recurrence, etc.)

8-SURGICAL TECHNIQUE

The placement of the spinal dynamic posterior stabilization device BDyn begins by the preparation of the pedicles of the vertebrae and by the pedicle aim and the fixation of the pedicle screws.

When tapping the vertebrae to prepare the hole before the insertion of the pedicle screw, flexion efforts on the tap are to be avoided. It must stay in the lined-up in order to make as straightest a pre-hole as possible to ensure as secure a screw placement as possible, with a straight guidance way. Thus, when tapping the vertebrae in order to stay in the line-up and to create a straight pre-hole.

Viewed in the sagittal plane, the pedicle screws must be as parallel as possible in order to let the biggest available space to place the cylindrical body of BDyn spinal shock absorber and facilitate its insertion between the head of the polyaxial pedicle screws. It is preferable to be sure that the adequate implant, with an appropriate size is used with the instrumentation which suit. A 2 mm space between the top of the cylindrical body of the BDyn spinal shock absorber and the lower part of the screw head must be respected to insure the shock absorption stock necessary for the good working of the BDyn device. For that, the surgeon can use the trial device RCBANTD50U to check if there is enough space to place the BDyn spinal shock absorber.

The mobile piston rod of the BDyn spinal shock absorber must be fixed in the head of the screw concentrically to the cylindrical body.

The surgery ends by the final tightening by using the dedicated torque limiting of the plugs on the head of the polyaxial pedicle screws by using systematically the suitable antirotator.

It is important to manipulate correctly the implants. The components of the implants to be arched do not have to be against-arched in the opposite direction. In any case they have to be chipped. These operations could effectuate the concentration of internal stresses which could become the place of an eventual failure of the implant.

9-PRESENTATION AND WARNING ABOUT ANCILLARIES

- The BDyn ancillaries are Class I medical devices, intended for temporary use and are re-usable.

Unlike the BDyn implants, the ancillaries included in the instrument set are furnished non sterile.

Before use, the ancillaries must be:

- cleaned by the appropriate mode of cleaning
- sterilized in an autoclave by moist heat

After use, the ancillaries have to be cleaned and decontaminated according to a suitable protocol especially by taking into account the reduction of the risk of transmission of unconventional transmissible agents - UTA.

In case of a patient suspected of having Creutzfeld-Jacob Disease (CJD), perform suitable decontamination.

In the event of a confirmed case of CJD, incinerate the ancillaries.

10-DECONTAMINATION AND CLEANING OF BDDYN ANCILLARIES

Decontamination and cleaning are performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments

The detachable instruments have to be dismantled

□Products advised for the cleaning

- an enzymatic detergent neutral
- Or products adapted for the cleaning (neutral or soft alkaline) of chirurgical instruments in compliance with the reglementation in effect.

- **Warning:** Don't use corrosive factor or caustic cleaning product

2 process are possible

A- Chart Automated cleaning in a disinfecter washer with manual Pre-cleaning:

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush (Do not utilize metal cleaning brushes)	Until all visible soils removed

Cleaning Steps	2	Pre-soak	Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.	<ul style="list-style-type: none"> - Detergent - Water jet pistol 	Until product is visually clean
	3	Ultrasonic cleaning	Clean the instrument in the ultrasonic tray including the solution. Be careful every unclench connector or screws and bolts by the vibrations will be reassure.	<ul style="list-style-type: none"> - Ultrasonic tray 	Minimum duration 10 minutes, this duration can vary according product
	4	Soak	Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position.	<ul style="list-style-type: none"> - Reserve osmosis water 	Minimum 30 seconds of soak
	5	Automated Washer	Place entire disassembled device into the automated washer	<ul style="list-style-type: none"> - Automated Washer - Disassembled Device 	Minimum total cycle time: 39 minutes

Instructions of cleaning in washer / disinfecter			
Cycles parameters	Time	Minimum temperature	Type of detergent/water
Pre-cleaning	2 minutes	<45°C / <113°F	Tap water
Cleaning	5 minutes	55°C / 131°F	Detergent
Draining	2 minutes	<45°C / <113°F	Tap water
Thermic draining	5 minutes	90°C / 194°F	Reverse osmosis water
Drying	25 minutes	> 70°C / > 158°F	Not applicable

Finishing Steps	6	Final Rinse	Rinse carefully with distilled water	<ul style="list-style-type: none"> - Distilled water 	Minimum 1 minute
	7	Final Drying	Dry devise utilizing medical quality filtered air	<ul style="list-style-type: none"> - Medical quality filtered air 	Until product is visually dry
	8	Tidying	The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization.		

B- Chart Just manual cleaning :

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	<ul style="list-style-type: none"> - Tap water - Soft bristle brush (Do not utilize metal cleaning brushes) 	Until all visible soils removed
	2	Drying	Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available	<ul style="list-style-type: none"> - Non-schedding wipe - Medical quality filtered compressed air 	Until product is visually dry
Disinfection Steps	3	Disinfection Application	Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument.	<ul style="list-style-type: none"> - Detergent 	Minimum duration 1 minute, this duration can vary according product
	4	Drying	After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water.	<ul style="list-style-type: none"> - Disposable cloth - Detergent 	Until product is visually dry
	5	Manual Disinfection	Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution.	<ul style="list-style-type: none"> - Soft bristle brush - Detergent - Syringe 	Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned
	6	Rinse and finition	Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water	<ul style="list-style-type: none"> - Reserve osmosis water - Syringe 	Minimum 1 minute
	7	Final Rinse	Rinse carefully with distilled water	<ul style="list-style-type: none"> - Distilled water 	Minimum 1 minute
	8	Final Drying	Dry devise utilizing medical quality filtered air	<ul style="list-style-type: none"> - Medical quality filtered air 	Until product is visually dry
	9	Tidying	The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization.		

Note: In case of **suspicion of Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers
In case of non respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

11-STERILIZATION OF BDYN ANCILLARIES

Sterilization is performed under the responsibility of the healthcare centers.
The following prescribed method and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).
In case of non respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for **which it will be fully responsible**.

1- The sterilization in autoclave is to make in a **specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

NB: It **isn't recommended** to realize the sterilization with peelables pouches **autoclavables** without a container Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

Steam sterilization parameters	Values
Cycle Type	Pre-vacuum
Set Point Temperature	134°C / 273°F
Cycle Time	Minimum 3 minutes
Dry Time	Minimum 30 minutes
Cool Time (inside and outside the autoclave chambers)	Minimum 40 minutes

3 - Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrality and the correct operation, of the locking systems (clipping, prehension) without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.

In case of non respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

12-STORAGE PRECAUTIONS

- The BDyn devices must be stored in a clean, dry and tempered place under atmospheric pressure, away from sunlight and rays in its original packaging. Ensure that the box is not dropped, hit or crushed.
- The ancillaries must be stored in a dedicated container or equivalent packaging to prevent any deterioration.

13-EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be done following recommendations of ISO 12891-1:2015 « Implants for surgery – Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be send back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package.

The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazards waste.

The elimination of a non-implanted device is not the subject of specific recommendations.

The surgeon will decide the necessity of the resection of the BDyn, its replacement or the rachis fusion according to the known and recognized techniques.

14-INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the healthcare center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time. For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim. Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

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Item	Patient information leaflet
1	<p>a. The name of the device : BDyn</p> <p>b. The model of the device : RCBDYMD55U, RCBDYSD55U</p>
2	<p>a. The intended purpose of the device : To restore the stabilization of the intervertebral thoracic, lumbar and sacral segments by preserving the anatomical lordosis and the deadening of the intervertebral joint.</p> <p>b. The kind of patient on whom the device is intended to be used : Adult human and female or male more than 40 kg.</p>
3	<p>Special operating instructions for the use of the device : The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique.</p>
4	<p>a. The intended performance of the device : The BDyn spinal shock absorber of the sterile spinal dynamic posterior stabilization device BDyn consists of a metallic cylindrical hollow part containing elastomer components made of silicone elastomers (polydimethyl siloxane (PDMS)) and long-term implantable polycarbonate urethane (PCU) which are bending out under the effect of a metallic piston rod connected with the vertebra of the treated segment by the pedicular screws tested and/or approved by the company COUSIN BIOTECH. The elastomer components assure the absorption of the mechanical loads applied on the intervertebral joint in compression, pulling, flexion-traction and lateral flexion. Some configurations enable fixation across several spinal segments.</p>

This leaflet was revised in March 2023

Page 59/64

Headquarters & Factory: COUSIN BIOTECH Allée des Roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –
Fax +33 (0) 3 20 14 40 13
www.cousin-surgery.com - **e-mail** f.pelletier@cousin-surgery.com

Item	Patient information leaflet
4	<p>b. Undesirable side effects : Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation: Infection, Pseudomeningocele, fistula, breach dura, persistent CSF leakage, meningitis, loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits, cauda equina syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss, urinary retention or loss of bladder control or other types of urological system compromise, scar formation possibly causing by a neurological compromise or compression around nerves and/or pain, fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery, herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery, canal adjacent stenosis, Non-union or pseudarthrosis, delayed union, mal union, Cessation of any potential growth of the operated portion of the spine, Loss of or increase in spinal mobility or function, Inability to perform the activities of daily living, Bone loss or decrease in bone density, Graft donor site complications including pain, fracture, or wound healing problems, Ileus, gastritis, bowel obstruction or loss of bowel control or other types of gastrointestinal system compromise, hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise, Reproduction system compromise, sterility, sexual dysfunction, Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc, Change in mental status, Death</p>
5	<p>When BDyn Should Not be Used : Active infectious process or significant risk of infection (immunocompromise), signs of local inflammation, Fever or leukocytosis, Morbid obesity, Pregnancy, Mental illness, Grossly distorted anatomy caused by congenital abnormalities, Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of the sedimentation rate unexplained by other diseases, elevation of the white blood count</p> <p>Suspected or documented metal allergy or intolerance, Any case where the implant components selected for use would be too large or too small to achieve a successful result, Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality, Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance, Any patient unwilling to follow postoperative instructions, Any case not describe in the indications, Traumas (i.e. fracture or dislocation), Abnormal curvatures (i.e. scoliosis and/or hyper lordosis), Tumors, Spondylolisthesis grade 2 and more, Pseudarthrosis and/or failed previous fusion, Severe bone resorption, osteomalacia, severe osteoporosis</p>

This leaflet was revised in March 2023

Page 60/648

Headquarters & Factory: COUSIN BIOTECH Allée des Roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –

Fax +33 (0) 3 20 14 40 13

www.cousin-surgery.com - **e-mail** f.pelletier@cousin-surgery.com

Item	Patient information leaflet
6	<p>MRI Information : The implants are composed of non-ferromagnetic materials and present a geometry non susceptible to generate induced currents. Moreover, as they are fixed to bone /or/ tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a safety use for MRI scan between 1.5 and 3 Tesla. As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential. It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams. For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.</p>
7	<p>Precautions : The BDyn devices are delivered sterile, sterilized by ethylene oxide. Before any use, inspect the integrity of the packaging and device (including peelable pouches). Do not use in the event of deterioration of the labels and/or the device and/or the packaging. Do not use if the device is out of date. The sterile spinal dynamic posterior stabilization device BDyn is only intended to adult human female or male more than 40 Kg. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. The installation of BDyn device can only be made with pedicular screws tested and/or approved by the company COUSIN BIOTECH. The eventual substitution of the BDyn spinal shock absorber has to be carried out only with the dedicated replacement rigid rods. You should never use implants made of stainless steel and titanium alloy in the same construct. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise the results. This device system is not intended to be the sole means of spinal support. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly and/or breakage of the device(s) will eventually occur. Preoperative and operating procedures, including precise knowledge of suitable surgical techniques, and proper selection of the good reference of the device adapted to the patient and its narrow setting up are important considerations in the successful use of the device by the surgeon. The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique Postoperative precautions</p>

This leaflet was revised in March 2023

Page 61/64


Headquarters & Factory: COUSIN BIOTECH Allée des Roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –
Fax +33 (0) 3 20 14 40 13
www.cousin-surgery.com - **e-mail** f.pelletier@cousin-surgery.com

Item	Patient information leaflet
8	<p>Materials included in the device : The BDyn Device components are manufactured from medical grade titanium TA6V (ISO 5832-3; ASTM F 136), silicone, long-term implantable polyurethane.</p> <p>BDyn implant: Titanium alloy Ti6Al4V ELI (ISO 5832-3; ASTM F 136), Polydimethyl siloxane (PDMS) MED 4770 unrestricted Polycarbonate urethane (PCU) Bionate® II 80A</p>
9	<p>Reporting a serious incident : If you encounter any serious incident with your BDyn implant you should report it immediately to your doctor. You should also report the incident to COUSIN BIOTECH via the link serviceclients@cousin-surgery.com or materiovigilance@cousin-biotech.com and to the Therapeutic Goods Administration using the following website: http://www.tga.gov</p>
10	<p>Legal Manufacturer & Manufacturing Facility : COUSIN BIOTECH, Allée des Roses, 59117 Wervicq-Sud, France</p>

This leaflet was revised in March 2023

Page 62/64

Headquarters & Factory: COUSIN BIOTECH Allée des Roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –
Fax +33 (0) 3 20 14 40 13
www.cousin-surgery.com - **e-mail** f.pelletier@cousin-surgery.com

en fr de it es	- Symbols used on labelling - Symboles utilisés sur l'étiquette - Bei Etiketten verwandete Symbole - Simboli utilizzati sull'etichetta - Símbolos utilizados en el etiquetaje	pt nl pl ru ελ cs	- Símbolos usados na etiqueta - Op de etikettering gebruikte symbolen - Symbole użyte na etykięcie - Символы, используемые на маркировке - Σύμβολα που χρησιμοποιούνται στις ετικέτες - Symboly použité na štítku	 www.cousin-biotech.com/ifu  IFU can be downloaded through the QR code. A hard copy can be sent within 7 days on request by email to ifurequest@cousin-biotech.com or by using the order form on our website.
----------------------------	---	----------------------------------	--	--

	LOT	REF		
en fr de it es pt nl pl ru ελ cs	Batch number Numéro de lot Chargennummer Numero di lotto Número de lote Número do lote Batchnummer Numer partii Номер партии Αριθμός παρτίδας Číslo šarže	Reference on the brochure Référence du catalogue Verweis auf der Broschüre Referimento del catalogo Referencia al foletto Referência do catálogo Referentie in de brochure Symbol katalogowy Ссылка на брошюру Αριθμός καταλόγου Odkaz v brožůře	Caution (See instructions for use) Attention, voir notice d'instructions Vorsicht (Siehe Gebrauchsanweisung) Attenzione (Vedere manuale istruzioni) Precaution (Veanse las instrucciones de uso) Atenção (Consultar nota de instruções) Vorsichtigheid (Zie gebruiksinstructies) Uwaga, patrz instrukcja obsługi Смотрите инструкцию по применению Δείτε τις οδηγίες χρήσης Viz pokyny k použití	Manufacturer Fabricant Hersteller Produttore Fabricante Fabricante Fabrikant Producent Производитель Κατασκευαστής Výrobce

	en fr de it es pt nl pl ru ελ cs	Do not re-use Ne pas réutiliser Nicht wiederverwenden Non riutilizzare No volver a utilizar Não reutilizar Niet hergebruiken Nie stosować ponownie Не предназначен для повторного использования Να μην επαναχρησιμοποιείται Nepoužívat znovu		en fr de it es pt nl pl ru ελ cs	Use before: year and month Utiliser jusque : année et mois Verbrauchen bis : Jahr und Monat Utilizzare entro e non oltre: anno e mese Úsese antes de: año y mes Utilizar até : ano e mês Te gebruiken vóór: jaar en maand Stosować do: rok i miesiąc Срок годности: год и месяц Ημερομηνία λήξης: έτος και μήνας Spotřebovat do: rok a měsíc
	en fr de it es pt nl pl ru ελ cs	Keep in a dry place Conserver au sec Trocken lagern Conservare in un luogo asciutto Manténgase en un lugar seco Manter em lugar seco Op een droge plaats bewaren Przechowywać w suchym miejscu Хранить в сухом месте Φυλάξτε το προϊόν μακριά από την υγρασία Uchovajte na suchém místě		en fr de it es pt nl pl ru ελ cs	Do not use if the packaging is damaged Ne pas utiliser si l'emballage est endommagé Nicht verwenden, wenn Packung beschädigt ist Non utilizzare se la confezione è danneggiata No usar si el embalaje está dañado Não usar se a embalagem estiver danificada Niet gebruiken indien de verpakking beschadigd is Nie używać, jeśli opakowanie jest uszkodzone Не использовать в случае повреждения упаковки Να μην χρησιμοποιείται, εάν η συσκευασία έχει υποστεί κάποια βλάβη Nepoužívejte je-li obal poškozený
	en fr de it es pt nl pl ru ελ cs	Do not re-sterilize Ne pas restériliser Nicht erneut sterilisieren Non ristilizzare No volver e esterilizar Não reesterilizar Niet opnieuw steriliseren Nie sterylizować ponownie Не подлежит повторной стерилизации Να μην επαναποστειρώνεται Znovu nesterilizujte		en fr de it es pt nl pl ru ελ cs	Keep away from sunlight A stocker à l'abri de la lumière du soleil Vor Sonnenlicht geschützt lagern Tenere lontano dalla luce Manténgase oscuras Manter ao abrigo da luz Vrij van licht houden Przechowywać w miejscu chronionym przed światłem Хранить в темном месте вдали от солнечных лучей Προφυλάξτε το προϊόν από το φως Uchovávejte mimo dosah slunečního záření

STERILE R	: Polyaxial pedicle screws, Replacements rigid rods	STERILE EO	: BDyn spinal shock absorber
en fr de it es pt nl pl ru ελ cs	Sterilized by gamma radiation. Produit stérile. Méthode de stérilisation : irradiation Sterilisierungsmethode : Gammastrahlung. Il prodotto è sterile. Sterilizzazione : raggi Gamma. Esterilizado por rayos gamma. Produto estéril. Método de esterilização: irradiação Sterilizováno gamma zářením. Steryliżowane przy użyciu promieni gamma. Стерилизовано гамма-излучением. Στείρο προϊόν. Μέθοδος αποστείρωσης: ακτινοβολία Sterilizované gamma žiarením	en fr de it es pt nl pl ru ελ cs	Sterilized by ethylene oxide Produit stérile. Méthode de stérilisation : oxyde d'éthylène Mit Ethylenoxid sterilisiert Sterilizzato all'ossido di etilene Esterilizado con óxido de etileno Produto estéril. Método de esterilização: óxido de etileno Gesteriliseerd met ethyleenoxide Steryliżowany tlenkiem etylenu Стерилизовавшийся с помощью окиси этилена Στείρο προϊόν. Μέθοδος αποστείρωσης: αιθυλενοξειδίο Sterilizované etylénoxidom