4D DOME[®]

SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

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CE
0297

NOT301_240215
Version du 15/02/2024

R_x only Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

Product under regulation 2017/745 (EU).

4DDOME®

SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT SINGLE USE STERILE PRODUCT

MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
ALLEE DES ROSES
59117 WERVICQ SUD
FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

The 4DDOME® device consists of a semi-resorbable dome and a semi-resorbable reinforcement mesh for anterior and posterior parietal reinforcement of the groin.

INDICATIONS

Repair and surgical reinforcement of inguinal and femoral hernias.

CONTRAINDICATIONS

Do not implant in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
- Infection
- Inflammation
- Fistula
- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes
- Recurrence
- Adhesion
- Erosion/Extrusion
- Seroma/lymphocele
- Hematoma
- Mesh migration
- Mesh deformation
- Allergic reaction
- Foreign body reaction
- Irritation nearby organ
- Obstruction

TARGET POPULATION

Adult male or female suffering from groin hernias with an indication of surgical repair with extraperitoneal implantation. These include:

- Inguinal hernias
- Femoral hernias

IMPLANTED MATERIALS

Polypropylene (PP) (non resorbable) - Poly-L-Lactic Acid (PLLA) (resorbable)

		Weight PP (g)	Weight PLLA (g)
4DDOM24GSR	Dome	0,03	0,26
	Flat mesh	0,21	0,64
4DDOM30GSR	Dome	0,05	0,41
	Flat mesh	0,21	0,64
4DDOM38GSR	Dome	0,08	0,60
	Flat mesh	0,21	0,64
4DDOME24SR	Dome	0,03	0,26
	Flat mesh	0,12	0,37
4DDOME30SR	Dome	0,05	0,41
	Flat mesh	0,12	0,37
4DDOME38SR	Dome	0,08	0,60
	Flat mesh	0,12	0,37

Origin neither human nor animal origin. Semi-resorbable

PERFORMANCES

The device 4DDOME® contributes to reduce symptoms from hernia and contributes to avoid complications of hernia.

LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:

https://www.cousin-biotech.com/d-6318526d8b749_under-construction.pdf

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is made with Polypropylene and Poly-L-Lactic Acid: non metallic components. Consequently, 4DDOME® is considered MR safe.

PRECAUTIONS FOR USE

The dome should always be installed with an anterior reinforcement.

The dome should not be cut.

4DDOME® devices must always be separated from abdominal cavity by peritoneum

4DDOME® devices are delivered sterile (ethylene oxide sterilisation).

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the mesh are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Before any use, inspect the integrity of packaging and device (of which blister / peel-away pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

Prevent the mesh from any contact with objects that may alter their surface.

This device must be implanted only by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As indicated on the product labelling, the 4DDOME® implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

STORAGE PRECAUTIONS

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

SURGICAL TECHNIQUE

4DDOME® Surgical Technique (1)

Indirect Hernia

① In indirect inguinal hernia, the sac is dissected up to the internal inguinal ring

② The 4DDOME® is placed in the defect maintaining the sac deeply. The 4DDOME® is then fixed by 1 or 4 stitches with non-resorbable suture to the resistant edge of the internal ring

③ Positioning the split onlay mesh around the spermatic cord and fixation with non-resorbable suture as in the Lichtenstein technique.

4DDOME® Surgical Technique (2)

Direct Hernia

① In direct hernia the sac is dissected up to the transversalis fascia.

② The 4DDOME® is placed in the opening whilst holding the sac deeply. The 4DDOME® is then fixed by 7 or 8 stitches with non-resorbable suture to the resistant edge of the transversalis fascia.

③ Positioning the split onlay mesh around the spermatic cord and fixation with non-resorbable suture as in Lichtenstein's technique.

EXPLANATION AND DISPOSAL OF DEVICES

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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