## SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

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Made in France

**NOT297\_ 240215**  
Version du 15/02/2024

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

**Product under regulation 2017/745 (EU).**

**SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT  
STERILE PRODUCT FOR SINGLE USE**

**MANUFACTURER IDENTIFICATION**

COUSIN BIOTECH  
ALLEE DES ROSES  
59117 WERVICQ SUD  
FRANCE

**IDENTIFICATION AND DESCRIPTION OF THE DEVICE**

4DFIX meshes are semi-resorbable parietal reinforcement implants.

**INDICATIONS**

4DFIX meshes are designed for the repair and reinforcement of inguinal and femoral hernias.

**CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
- Recurrence
- Hematoma
- Foreign body reaction
- Infection
- Adhesion
- Mesh migration
- Irritation nearby organ
- Inflammation
- Erosion/Extrusion
- Mesh deformation
- Obstruction
- Fistula
- Seroma/lymphocele
- Allergic reaction
- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic)

**TARGET POPULATION**

Adult male or female suffering from groin hernias with an indication of surgical repair with extraperitoneal implantation. These include:

- Inguinal hernias
- Femoral hernias

**IMPLANTED MATERIALS**

Polypropylene (PP) (non resorbable) – Poly L Lactic Acid (PLLA) (resorbable). Blue thread (depending on references): monofilament polypropylene (non resorbable), dye: [phthalocyaninato(2-)] copper.

	Weight of PLLA (g)	Weight of PP (g)	Weight of the blue colorant [Phthalocyaninato(2-)] copper (mg)
<b>VCBFXPRSRU</b>	0,77	0,27	< 0,01
<b>VCBFXPRSLU</b>	0,77	0,27	< 0,01
<b>VCBFXPRMRU</b>	0,99	0,35	< 0,01
<b>VCBFXPRMLU</b>	0,99	0,35	< 0,01
<b>VCBFXPRLRU</b>	1,09	0,38	< 0,01
<b>VCBFXPRLLU</b>	1,09	0,38	< 0,01
<b>VCBFXRABAU</b>	1,01	0,36	NA (no landmark)

VCBFX1317U	1,21	0,42	< 0,01
VCBFX1215U	1,05	0,37	< 0,01
VCBFXLAPLU	0,98	0,34	< 0,01
VCBFXLAPRU	0,98	0,34	< 0,01
VCBFX1717U	1,97	0,69	NA (no landmark)
VCBFX812FU	0,51	0,18	NA (no landmark)
VCBFX613LU	0,47	0,17	NA (no landmark)
VCBFX613RU	0,47	0,17	NA (no landmark)
VCBFX0812U	0,64	0,23	NA (no landmark)

No human nor animal origin. Partially absorbable.

### **PERFORMANCES**

The device 4DFIX contributes to reduce symptoms from hernia and contributes to avoid complications of hernia.

### **LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:**

[https://www.cousin-biotech.com/d-6318526d8b749\\_under-construction.pdf](https://www.cousin-biotech.com/d-6318526d8b749_under-construction.pdf)

### **MRI SAFETY**

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with Polypropylene and Poly-L-Lactic Acid: non-metallic components. Consequently, 4DFIX is considered MR safe.

### **PRECAUTIONS FOR USE**

4DFIX must always be separated from abdominal cavity by peritoneum. The rough side with loops should be placed facing the inguinal floor (see the surgical technique below). ~~abdominal wall.~~

4DFIX prosthesis are delivered sterile (ethylene oxide sterilisation).

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the mesh are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

Prevent the mesh from any contact with objects that may alter their surface.

Meshes for laparoscopic surgery have a blue landmark that should be positioned towards the medial side of the patient (see the surgical technique below).

This device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

### **IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As specified on the product's labelling, the 4DFIX mesh is for single use only. It cannot be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence).

### **INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON**

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

### **STORAGE PRECAUTIONS**

To be stored in a dry place away from light and at room temperature in its original packaging.

### **SURGICAL TECHNIQUE**

#### **DESCRIPTION OF THE SURGICAL PROCEDURES:**

##### **Open Inguinal hernia repair : Lichtenstein**

- The patient is placed in the supine position
- Local or general anaesthesia
- Iliac incision of approximately 5 cm
- Dissection is performed accordingly to the Lichtenstein technique
- The hernia sac is identified and reduced
- The transversalis fascia can be tightened with a non-absorbable suture. With direct hernias it helps to avoid any immediate post-operative recurrence

- 4DFIX mesh can be cut to fit the patient's anatomy
- The prosthesis is progressively put in place
  - Begin by placing 4DFIX on the pubic bone and with the rough side down (facing the inguinal floor)
  - The upper part of the prosthesis is placed on the conjoint tendon
  - Carry placing the prosthesis towards the deep inguinal ring so that the first arm is placed behind
  - Then the second arm is placed so that it encircles the spermatic cord
  - Fixation of the 4DFIX mesh with non-absorbable sutures, according to surgeon's technique
  - The spermatic cord rests on the prosthesis
  - Suction drain can be placed beneath the external oblique aponeurosis, especially in large inguinal hernias, where an extensive dissection was performed
  - The aponeurosis of external oblique is then closed with absorbable sutures
  - Before closure of the surgical incision, its edges are infiltrated with a long-acting local anaesthetic, such as Naropein®
  - Prophylactic antibiotics is usually given preoperatively

### **Laparoscopic procedure (TAPP or TEP)**

#### 1/ Example of a TAPP repair

- General anaesthesia
- A pneumoperitoneum is created in the usual fashion
- Intra-abdominal findings are reported [intra-abdominal pathology and inguinal hernia orifice and sac].
- The two additional trocars are inserted under direct vision.
- Dissection is initiated and carried out following classical TAPP technique
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected.
- The 4DFIX mesh is inserted through the Trocar into the intra-abdominal cavity and deployed with the rough side down, on the inguinal region. The blue thread is the anatomic landmark for the pubic bone. The 4DFIX mesh can be attached to Cooper's Ligament, around and lateral to the Inferior Epigastric Vessels using glue, tacks, staples, or other mesh fixating device, according to Surgeon's usual technique. Decision regarding fixating the mesh or not, remains the surgeon's decision. The operator should be meticulous to avoid the triangle of doom, triangle of pain, iliac vessels and to place fixations lateral to the inguinal ring. Mesh fixation points should be at least 1cm from the edge of the mesh with 1 cm spacing between fixation points.
- The operator should check the 4DFIX mesh is well anchored to the surrounding structures.
- The peritoneum is closed meticulously and no defect between the peritoneum and the abdominal wall should be left open. In addition, it should cover the entire Mesh.
- Trocars are removed under direct vision. Fascia of the sub-umbilical trocar site is closed as needed.

#### 2/ Example of a TEP repair

- General anaesthesia
- Create the pro-peritoneal space.
- The 10 mm trocar is inserted.
- Create the pneumo-pro-peritoneum.
- Two 5 mm ports are inserted under direct vision.
- If a peritoneal defect is created at this time it should be immediately closed.
- The anatomical landmarks are identical to when performing a TAPP repair.
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected. This is done by gentle traction with atraumatic graspers.
- The 4DFIX mesh is inserted into the pre-peritoneal space and deployed over the inguinal region. The blue thread is the anatomic landmark for the pubic bone.
- 4DFIX mesh is first secured to Cooper's Ligament, and the superior aspect of the pubic ramus. The 4DFIX mesh is then anchored on the posterior aspect of the Linea Alba. This anchoring is continued around and lateral to the Inferior Epigastric Vessels. A few fixations may be used to fix 4DFIX lateral to the internal ring. Mesh fixation or not, remains the surgeon's decision.
- The 4DFIX mesh should gently and generously cover the Iliac Vessels without major gaps.
- The operator should check the Mesh is well anchored to the surrounding structures. The Mesh should not migrate and remain in place.
- The pre-peritoneal space will be checked for any peritoneal defects. If any, they should be closed. The 5 mm trocars are removed under direct vision. The pre-peritoneal space will collapse.
- If there are any questions about a missed peritoneal defect, a completion laparoscopy should be performed.
- Trocars are removed and the fascial defect closed with the appropriate suture. The skin edges are approximated in the usual manner.

### **Size of the prosthesis:**

The 4DFIX product range is composed of meshes of various sizes and shapes:

- Pre-cut meshes with keyhole for LICHTENSTEIN inguinal and hernia repair
- Oval, pre-shaped and rectangular meshes for TAPP/TEP inguinal hernia repair

### **EXPLANATION AND DISPOSAL OF DEVICES**

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following

the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste. There are no specific recommendations for the disposal of a non-implanted device.

### **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

### **INFORMATION REQUEST AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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